

IYASOCARE®



FOCUS
RESEARCH
DEVELOPMENT
PRODUCTION

📍 Chongqing Iyasocare Medical Co.,Ltd.
📍 Free Trade Tower, Xiyong street. Chongqing, China
✉️ iyaso@iyasocare.com
☎️ 0086-136 1835 0425


🌐 www.iyasocare.com
iyasocare.en.alibaba.com

CHONGQING IYASOCARE MEDICAL CO.,LTD.

PRODUCT CATALOG BROCHURE

COMPANY INTRODUCTION

Chongqing IYASOCARE MEDICAL CO.,LTD. was established in 2006 and is headquartered in the beautiful mountain city of Chongqing in China. It is an integrated industry and trade company specializing in medical device R&D, production, sales and consulting services, with business all over the world. The company's main products include wheelchairs, hospital beds, electric scooters, walkers, commode chair, etc.; the company is committed to providing customers with reasonable, complete, reliable, and high-quality medical procurement programs; since the company was founded, it has served the United States, the United Kingdom, and Australia, France, Saudi Arabia, Egypt, Pakistan, Japan, South Korea and other global customers; and has obtained medical device registration certificate, medical device production license, medical device business license, ISO13485, CE, TUV, MDR, CFS and other qualifications Certificate; is a supplier of continuous strategic cooperation with China Medical Device Association.



AUTHENTICATION CERTIFICATE

Medical Management System Certification
Award
CHONGQING IYASOCARE MEDICAL CO.,LTD.
Room 305-306, Building 8, Chongqing Railway Post Office Logistics and Storage Project,
Xiyang Group 15th Floor, Xiyang Street, Shapingba District, Chongqing

The Medical management system of the above organizations has been reviewed and complied with **ISO 13485:2016** Requirements of applicable clauses of Medical management system standards

This system is valid for the **Production and sales of Medical equipment, Wheelchair, Manual wheelchair, Electric wheelchair, Scooter, Electric handcycle, Manual bed, Electric Bed, Medical bed products**

Under the condition that the certificate holder's management system continues to meet the management system standards, the certification deadline is three years. And query on www.ccc-iso.com. The validity of this certificate must be confirmed by the center through regular supervision and review.

Issue date: December 26, 2023 Valid period: November 25, 2028

ISIRI
IAF
UKAS
France
Spain European Certified organization Limited
www.ccc-iso.com

CE DECLARATION OF CONFORMITY
ACCORDING TO (EU) 2017/745 MEDICAL DEVICE REGULATION

Manufacturer
Name: Chongqing Iyascare Medical Co., Ltd
Address: 406, 4/F, Software Park Phase 2, 8 Xiyuan North Street (Free Trade Tower),
The Netherlands
SRN: NL-AR-00000247

Product Information
Name: ELECTRIC WHEELCHAIR
Model: see annex
EMDN: Y122127
Basic UDI-DI: /
Classification: Class I, According to Rule 13, Annex VII, Regulation (EU) 2017/745

Declaration
We herewith declare that the above-mentioned products meet the requirements of Medical Device Regulation (EU) 2017/745 and the applicable standards above.

Signature: [Signature] Date: 2023.12.26
Position: GM Place: Chongqing/China

CE DECLARATION OF CONFORMITY
ACCORDING TO (EU) 2017/745 MEDICAL DEVICE REGULATION

Manufacturer
Name: Chongqing Iyascare Medical Co., Ltd
Address: 406, 4/F, Software Park Phase 2, 8 Xiyuan North Street (Free Trade Tower),
The Netherlands
SRN: NL-AR-00000247

Product Information
Name: Electric scooter
Model: see annex
EMDN: Y122124
Basic UDI: /
Classification: Class I, According to Rule 13, Annex VII, Regulation (EU) 2017/745

Declaration
We herewith declare that the above-mentioned products meet the requirements of Medical Device Regulation (EU) 2017/745 and the applicable standards above.

Signature: [Signature] Date: 2023.12.28
Position: GM Place: Chongqing/China

Certificate of Conformity

Certificate No.: ECT2021P095388R-1
Applicant: CHONGQING IYASOCARE MEDICAL CO.,LTD.
Room 201-205,115 Longan Avenue, Tuzhu Town, Shapingba District, Chongqing Province, China
Product Name: Wheelchair, Manual wheelchair, Electric wheelchair, Scooter, Electric handcycle, Commode chair, Walking aids, Transfer wheelchair

Trade Name: N/A
Model: (see the following annex)
Classification: Class I (according to the Manufacturer's declaration)
Review goal: Verification of the presence of Technical Documentation compatible with the Medical Devices Regulation (EU) 2017/745

This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the Technical Documentation shared with us by the manufacturer is compatible with the European Standard for Medical Devices Regulation.

The manufacturer is responsible for CE Marking process and not exempted to carry out all necessary compliance activities. This document has been issued on the basis of the regulation on ECT Voluntary Mark for the certification of products.

Shenzhen ECT Testing Technology Co., Ltd.
3th Floor, EverBuilding Business Center, Building A, No. 26 Gonghe Industrial Road, Shenzhen, Shenzhen City, China.
Tel: +86-755-87587803 Fax: +86-755-87587803 <http://www.ectlab.com>

TÜVRheinland Made-in-China

Report Number: MIC-AS129874 Audit Date: 27 Aug. 2022

This report is issued by Focus Technology Co., Ltd (Made-in-China.com) and the supervising inspectorate (TUV Rheinland) to confirm that:

Company Name: Chongqing Iyascare Medical Co., Ltd.
Showeroom: <https://iyascare.com/made-in-china>
Address: Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China
Product: Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products

has been on site audited for the following scope of activity:

- General Information
- Foreign Trade Capacity
- Product Research & Development Capacity
- Management System and Product Certification
- Production Capacity & Quality Control
- The actual content of the part relates only to its associated enterprise & Working Environment
- Energy Saving and Emission Reduction
- Photo

General Comments:
Chongqing Iyascare Medical Co., Ltd. was founded in 2010, is a trading company with 10 employees, the company is located at Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China. The company mainly supply Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products. The company has successful foreign trading experience in North America, South America etc.

Signature: [Signature] Date: 16 Feb 2024
On behalf of TÜVRheinland

CEIG Ministerie van Volksgezondheid, Welzijn en Sport

Report Number: MIC-AS129874 Audit Date: 27 Aug. 2022

This report is issued by Focus Technology Co., Ltd (Made-in-China.com) and the supervising inspectorate (TUV Rheinland) to confirm that:

Company Name: Chongqing Iyascare Medical Co., Ltd.
Showeroom: <https://iyascare.com/made-in-china>
Address: Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China
Product: Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products

has been on site audited for the following scope of activity:

- General Information
- Foreign Trade Capacity
- Product Research & Development Capacity
- Management System and Product Certification
- Production Capacity & Quality Control
- The actual content of the part relates only to its associated enterprise & Working Environment
- Energy Saving and Emission Reduction
- Photo

General Comments:
Chongqing Iyascare Medical Co., Ltd. was founded in 2010, is a trading company with 10 employees, the company is located at Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China. The company mainly supply Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products. The company has successful foreign trading experience in North America, South America etc.

Signature: [Signature] Date: 16 Feb 2024
On behalf of CEIG

Electric Wheelchair

>>> YSE123

Classic model



FEATURES:

1. Lightweight folding design, easy to carry and travel
2. Switch to manual model(3s) or manual switch
3. Detachable footrest
4. Drop back handle
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lead Acid
 Cushion: Black/Red/Blue

>>> YSE111



>>> YSE113



>>> YSE112



>>> YSE122

FDA

SFDA



>>> YSE126



>>> YSE131



MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
44.5KG/49KG	100KG	118CM	63CM	93CM	45CM	45CM	10 inches	16 inches	39CM	50CM	72*38*74CM	91*45*84CM

Electric Wheelchair

YSE109

Manual reclining



FEATURES:

1. Switch to manual model
2. Easy to fold
3. Reclining backrest 90°-150°
4. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15 km
 Max Speed: 6km/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lead Acid
 Parking Brake: Electronic parking brake
 Cushion: Red/Black/Blue

YSE108

Electric reclining



YSE118

FDA

SFDA

Electric reclining



MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
50KG/56KG	130KG	110CM	64CM	125CM	46CM	45CM	10 inches	16 inches	75CM	54CM	79*39*74CM	91*45*84CM

YSE130M

Aluminum



FEATURES:

1. It can be carried on a plane
2. Controller charging
3. Lightweight folding design, easy to carry and travel, easy manual folding
4. Switch to manual model
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15-25KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lithium battery
 Parking Brake: electromagnetic brake
 Cushion: Black

YSE130X

Auto folding



YSE130XB

Electric reclining



MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
29KG/ 33KG	120KG	98CM	63CM	110CM	48CM	45CM	8 inches	16 inches	43CM	50CM	78*48*63CM	67*50*80CM

Electric Wheelchair

YSE110

Heated seats

Ventilated seats



FEATURES:

1. Seat and backrest can be massaged, heated and ventilated, three-level adjustment
2. Headrest with Bluetooth audio connection
3. Aluminum alloy lightweight folding frame
4. Foldable armrest and footrest.

PARAMETERS:

Charge time: 4-6H
 Driving Range:15KM
 Max Speed:4.5KM/h
 Climbing ability:12°
 Motor:180*2 brushless motor
 Battery:24V 11Ah Lithium
 Parking Brake: electromagnetic braking
 Cushion:Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
25/28KG	100KG	100CM	62CM	96CM	45CM	45CM	8 inches	8 inches	44CM	50CM	67*41*62CM	65.5*44*74CM

YSE121

FDA

SFDA



FEATURES:

Charge time: 4-6H
 Driving Range: 10-20KM
 Max Speed: 6KM/h
 Climbing ability:12°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lithium battery
 Parking Brake:electromagnetic brake
 Cushion: Black

PARAMETERS:

- 1.It can be carried on a plane
- 2.Controller charging
- 3.Lightweight folding design, easy to carry and travel ,easy manual folding
- 4.Switch to manual model
- 5.Flip-up armrest

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
29KG/33KG	120KG	95CM	63CM	98 CM	45CM	45CM	8 inches	12 inches	43CM	55CM	78*63*48CM	87*66*44 CM

YSE120

Aluminum

Carbon fiber



FEATURES:

Charge time: 4-6H
 Driving Range: 10-20KM
 Max Speed: 6KM/h
 Climbing ability:9°
 Motor: 24V 180W*2Pcs brushless motor
 Battery: 24V 12Ah lithium battery
 Parking Brake:electromagnetic brake
 Cushion: Black

PARAMETERS:

- 1.It can be carried on a plane
- 2.Controller charging
- 3.Lightweight folding design, easy to carry and travel ,easy manual folding
- 4.Switch to manual model
- 5.Flip-up armrest

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
20K/22KG	100KG	95CM	63CM	98 CM	44CM	45CM	7 inches	11 inches	/	45CM	78*63*48CM	63*40*83.5CM

YSE124

Aluminum



FEATURES:

- 1.Lightweight folding design,easy to carry and travel
- 2.Easy to portable
- 3.Detachable footrest
- 4.Drop back handle
- 5.manual or electric model switch

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15-25KM
 Max Speed: 6KM/h
 Climbing ability:15°
 Motor: 250W*2Pcs/ 150W*2 brushless motor
 Battery: 24V 12Ah lithium battery
 Parking Brake: electric parking brake
 Cushion:Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
20KG/24KG	100KG	95CM	58CM	89CM	43CM	50CM	7 inches	10 inches	39CM	50CM	72*38*74CM	62*35*75CM

Electric Wheelchair

YSE125

Steel



FEATURES:

1. Lightweight folding design, easy to carry and travel
2. Easy to portable
3. Detachable footrest
4. Drop back handle
5. manual or electric model switch

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15-25KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lead batteries
 Parking Brake: electric parking brake
 Cushion: Black/Red/Blue

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	95CM	58CM	89CM	43CM	45CM	7 inches	10 inches	50CM	50CM	60*29*74CM	62*35*75CM

YSE130

Steel



FEATURES:

1. It can be carried on a plane
2. Controller with USB mobile phone charging
3. Lightweight folding design, easy to carry and travel ,easy manual folding
4. Switch to manual model
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15 KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lead Acid
 Parking Brake: electromagnetic brake
 Cushion: Black(8CM thickened)

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	94CM	62CM	89CM	44CM	45CM	8 inches	12 inches	43CM	50CM	75*62*34CM	70*43*79CM

YSE128

Steel



FEATURES:

1. Lightweight folding design, easy to carry and travel
2. Switch to manual model
3. Detachable footrest
4. Drop back handle
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15 KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 12V 12Ah (6*2 lead Acid)
 Parking Brake: electromagnetic brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	92CM	63CM	85CM	46CM	43CM	8 inches	12 inches	50CM	50CM	70*36*71CM	71*43*76CM

YSE130L

Aluminum



FEATURES:

1. It can be carried on a plane
2. Controller with USB mobile phone charging
3. Lightweight folding design, easy to carry and travel ,easy manual folding
4. Switch to manual model
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15-25KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lithium battery
 Parking Brake: electromagnetic brake
 Cushion: Black(8CM thickened)

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	94CM	62CM	89CM	44CM	45CM	8 inches	12 inches	43CM	50CM	75*62*34CM	70*43*79CM

Electric Wheelchair

YSE134



FEATURES:

1. Universal controller
2. Electric strong power offroad wheels wheelchair for disabled and elderly people.
3. Anti-slip design, the wheelchair will not slide backwards on the slope
4. Modular design, the whole disassembly, can put the auto trunk
5. Safe low-speed start design to avoid potential safety hazards caused by excessive acceleration
6. The seat can go up 45 degrees, and the back can lie down 160 degrees with electric adjustment.
7. The foot rest can go up 70 degrees adjustable.
8. Seat rotation design, 360° seat free adjustment.

PARAMETERS:

Charge time: 8-10H
 Driving Range: 20-40 KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Battery: 24V 15Ah/30Ah lithium
 Brake type: Electromagnetic brake
 Braking distance ≤800mm
 Obstacle crossing height ≤80mm
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	115CM	58CM	100CM	45CM	42CM	8 inches	12 inches	/	44CM	/	/

YSE144

Ergonomic seat



FEATURES:

1. It can be carried on a plane
2. Controller with USB mobile phone charging
3. Lightweight folding design, easy to carry and travel, easy folding. Switch to manual model
4. Flip-up armrest
5. Ergonomic chair

PARAMETERS:

Charge time: 8-10H
 Driving Range: 20 KM
 Max Speed: 6KM/h
 Climbing ability: 8°
 Motor: 250W*2Pcs brushless motor
 Battery: 24V 12Ah lithium
 Brake type: Electromagnetic brake
 Braking distance ≤700mm
 Obstacle crossing height ≤40mm
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
31.7KG/39KG	120KG	100CM	65CM	105CM	44CM	45CM	7 inches	12 inches	/	47-51CM	65*36*77CM	69*40*86CM

YSE149



FEATURES:

1. 10" McNamee wheels 13" inflatable rubber rear wheel
2. Flying wing shock absorption system
3. Automotive grade seat cushion, ergonomic design, high-end leather
4. Intelligent programmed brush controller
5. Fully automatic electromagnetic brake

PARAMETERS:

Charge time: 4-6H
 Driving Range: 15-25km
 Max Speed: 6km/h
 Climbing ability: 10°
 Battery: 24V 12/ 20A Lithium battery
 Motor: 250W x2
 Parking Brake: Electromagnetic brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
50/59KG	130KG	94.5CM	64CM	88CM	47CM	41CM	10 inches	13 inches	/	55CM	94.5*64*66CM	97*67*45 CM

YSE150



FEATURES:

Seat width: 40 cm
 Seat depth: 44cm
 Seat height: 54 cm
 Product size: 102*55*99 cm
 Rear wheel size: 9 inches
 Front wheel size: 8 inches
 N/G: 30/ 35KG
 Frame material: steel
 Load: 130 KG
 Packing size:

PARAMETERS:

1. 8-inch McNamee front wheel
2. aluminum alloy frame
3. Seat height adjustable, telescopic footrest, foldable armrests
4. Intelligent controller with lighting function
5. Fully automatic electromagnetic brake
6. 3D butterfly wing cushion backrest

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
30/35KG	130KG	102CM	55CM	99CM	40CM	44CM	8 inches	9 inches	/	54CM	/	/

Electric Wheelchair

YSE154

Steel



FEATURES:

- 1, Folding front and back
- 2, Armrests can be turned over
- 3, Pedal can be turned over
- 4, Battery quick disassembly
- 5, Controller front and back adjustment, left and right interchangeable

PARAMETERS:

Charge time: 4-6H
 Driving Range: 15KM
 Max Speed: 7KM/h
 Climbing ability:12°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lead Acid
 Parking Brake: manual brake and electric parking brake
 Cushion: Black/Red/Blue

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
33/36kg	100KG	102CM	67CM	89M	43CM	44CM	8 inches	12 inches	78CM	51CM	70*42*73CM	88*42*86CM

IYASOCARE®

Off Road Wheelchair

YSE101

All terrain



FEATURES:

1. Universal controller.
2. Electric strong power offroad wheels wheelchair for disabled and elderly people.
3. Can climb 10cm high stairs.
4. Anti-slip design, the wheelchair will not slide backwards on the slope.
5. Standard USB charging port, safety belt.

PARAMETERS:

Charge time: 8-10H
 Driving Range: 30-60km
 Max Speed: 0-8 km/h
 Climbing ability: 15°
 Motor: 24V 320W*2
 Battery: 24V 40A lead acid
 Cushion: Red/Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
65/63KG	200KG	115CM	66CM	110CM	45CM	45CM	10 inches	14 inches	78CM	56CM	/	115*65*60CM

YSE102

All terrain



FEATURES:

1. Electric strong power offroad wheels wheelchair for disabled and elderly people.
2. Can climb 10cm high stairs.
3. Anti-slip design, the wheelchair will not slide backwards on the slope.
4. Standard USB charging port, safety belt
5. Speed and battery life will vary with the battery & motor power changes.
6. The seat can be adjusted 170 degrees to lie flat.

PARAMETERS:

Charging time: 8-10h
 Driving range : 35-75km
 Top speed : 20 km/h
 Climbing capacity : 15°

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
80/92 KG	200KG	129CM	70CM	124 CM	45CM	45CM	13 inches	14 inches	78CM	56CM	/	115*65*60CM

YSE101A

All terrain



FEATURES:

1. Upgrade the thickened and widened seat, made of top grade leather.
2. Upgrade the new thickened leather armrests.
3. The body volume is narrower and shorter.
4. The rear wheel hub has increased the night driving lights and turn signals, more secure.
5. Foot rest can be adjusted forward and backward.

PARAMETERS:

Charge time: 6-8H
 Driving Range: 25-60km
 Max Speed: 10 km/h
 Climbing ability: 15°
 Motor: 24V 250W*2
 Battery: 24V 40Ah lead acid

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
80/85KG	200KG	110CM	62.5CM	121CM	45.5CM	45CM	10 inches	14 inches	/	56CM	/	115*65*60CM

YSE102A

All terrain



FEATURES:

1. Upgraded 800w motor, max speed 25km/hour.
2. Upgraded disc brakes on both front and rear, and added oil brakes on the rear wheel.
3. Upgraded wider and thicker leather seats than the first generation.
4. Upgraded the front and rear wheels.
5. Upgraded the center integrated pedal.

PARAMETERS:

Charging time: 6-8h
 Driving range: 25-75km
 Max speed: 25 km/h
 Climbing capacity: 15°
 Motor: 800W differential motor
 Battery: 48v 20ah / 32ah lead acid / 60v 30ah lithium battery
 Parking brake: disc brake, oil brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
105/110KG	200KG	126CM	65.5CM	125CM	46CM	45CM	14 inches	16 inches	/	57CM	/	115*65*64CM

Off Road Wheelchair

YSE103

All terrain

Double seat



FEATURES:

1. Increase the rear seat, can be traveled by two people.
2. Dual handles are easy to control, the tires are enlarged and widened, equipped with spring shock absorption.
3. Large motor power, fast speed, large battery capacity, long range.
4. Suitable for paraplegic, hemiplegic patients, lower limb disabled and mobility impaired elderly.
5. Suitable for long-distance traveling, demanding speed and dealing with poor road surface.

PARAMETERS:

Charging time: 6-8h
 Driving range :25-75km
 Max speed: 20 km/h
 Climbing capacity :12°
 Motor: 800W differential motor
 Battery: 48v 20ah / 32ah lead acid / 60v 30ah lithium battery
 Parking brake: drum brake
 Cushion:Red/Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	200KG	166CM	90CM	125CM	46CM	45CM	13 inches	14 inches	/	57CM	/	145*80*70CM

YSE105

All terrain



FEATURES:

1. Electric strong power offroad wheels wheelchair for disabled and elderly people.
2. Standard USB charging port, safety belt, optional selfie stick bracket and personalized storage basket.
3. Speed and battery life will vary with the battery & motor power changes.
4. Four driving modes are available: off-road, climb, normal, and indoor.
5. 180° reclining adjustable car seat.

PARAMETERS:

Charging time: 4-6h
 Driving range :40-60km
 Top speed :28 km/h
 Battery: 48v 20ah lead acid
 Parking brake: disc brake

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	200KG	100CM	63CM	58CM	49CM	46CM	14 inches	14inches	/	58CM	/	135*70*90CM

YSE104

All terrain

Electric reclining



FEATURES:

1. Universal controller
2. Electric strong power offroad wheels wheelchair for disabled and elderly people.
3. Can climb 10cm high stairs.
4. Anti-slip design, the wheelchair will not slide backwards on the slope
5. Standard USB charging port, safety belt.

PARAMETERS:

Product size:115*68*110 cm
 Rear wheel size: 14 inches
 Front wheel size: 10 inches
 Frame material: steel
 Load:200KG
 Packing size:115*65*60cm
 Wooden box:120*72*79cm 140KG

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	200KG	115CM	68CM	110CM	45 CM	45CM	10 inches	14 inches	78CM	52CM	/	115*65*60CM

YSE107

All terrain



FEATURES:

1. Suitable for disabled and elderly people's electric powerful cross-country wheelchair.
2. It can climb stairs of 10 centimeters high.
3. Anti-slip slope design, the wheelchair will not slide backward or roll over on the slope.
4. Comes standard with USB charging port, leather seat, seat belt, turn signal, headlight and horn.
5. Speed and battery life vary with battery and motor power.
6. Three driving modes are available: off-road, normal, and indoor.

PARAMETERS:

Charging time: 8-10h
 Driving range :40-60km
 Top speed :28 km/h
 Climbing capacity :15°
 Battery: 48v 20ah lead acid
 Parking brake: electromagnetic brake
 Cushion:black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	200KG	125CM	64CM	108 CM	49CM	46CM	14 inches	14 inches	/	58CM	/	135*70*90CM

Small Light Electric Wheelchair

YSE106

FDA



FEATURES:

- 1.It can be carried on a plane
- 2.18kg Lightweight folding design, easy to carry and travel
- 3.Intelligent electronic brake
- 4.Flip-up armrest
- 5.360° Adjustable push handle

PARAMETERS:

Charge time: 6-8H
 Driving Range: 10-15KM
 Max Speed:8KM/h
 Climbing ability:6°
 Motor:250W*2 Double hub motor
 Battery:24V 6Ah Lithium
 Parking Brake>manual brake and electronic parking brake
 Cushion:Black/blue

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	93CM	58CM	90CM	36CM	36CM	6 inches	8 inches	43CM	53CM	80*63*40CM	68*58*35CM

YSE129

Aluminum



FEATURES:

1. High-strength aluminum alloy frame, sturdy and durable
2. Brushless electromagnetic brake motor, safe and no slippery slope
3. Low noise Li-ion battery, lightweight and convenient
4. Long life brushless motor and controller, 360 degree flexible control.
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range:15 KM
 Max Speed:6KM/h
 Climbing ability:13°
 Motor:150W*2Pcs brushless motor
 Battery:24V 6Ah/12Ah lithium battery
 Parking Brake:electric brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	96CM	59CM	90CM	/	44CM	7 inches	10 inches	43CM	50CM	69*59*32CM	71*34*61CM

YSE119

Aluminum



FEATURES:

1. High-strength aluminum alloy frame, sturdy and durable
2. Brushless electromagnetic brake motor, safe and no slippery slope
3. Low noise Li-ion battery, lightweight and convenient
4. Long life brushless motor and controller, 360 degree flexible control.
5. Flip-up armrest

PARAMETERS:

Charge time: 4-6H
 Driving Range:15 KM
 Max Speed:6KM/h
 Climbing ability:8°
 Motor: 180w*2 brushless motor
 Battery:24V 6Ah Lithium
 Parking Brake>manual brake and electronic parking brake
 Cushion:Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	112CM	60CM	94CM	43CM	45CM	7 inches	12 inches	44CM	48CM	40*60*84CM	43*62*87CM

YSE132

Shock absorbers



FEATURES:

- Charge time: 4-6H
 Driving Range:5-10KM
 Max Speed:6KM/h
 Climbing ability:9°
 Motor:250W*2Pcs
 Battery:24V 6Ah/ 10Ah/ 20Ah lithium battery
 Parking Brake:electromagnetic brake
 Cushion: Black

PARAMETERS:

- 1.Aluminum alloy frame.
- 2.Breathable backrest, high quality leather padding.
- 3.Front wheel shock absorber
- 4.Rear flip armrest is equipped with PU armrest pad.
5. Front and rear wheels are equipped with PU non-inflatable honeycomb wheels.
- 6.Adjustable leg height

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	117CM	53CM	93M	41CM	40CM	8 inches	10 inches	/	48CM	/	55*34*82 CM

Small Light Electric Wheelchair

YSE136

Ergonomic seat



FEATURES:

Charge time: 4-6H
 Driving Range: 5-10KM
 Max Speed: 6KM/h
 Climbing ability: 9°
 Motor: 250W*2Pcs
 Battery: 24V 6Ah/ 10Ah/ 20Ah lithium battery
 Parking Brake: electromagnetic brake
 Cushion: Black

PARAMETERS:

1. Aluminum alloy frame, sand black finish.
2. Ergonomic seat plate with removable seat cushion.
3. Front wheel shock absorber.
4. Rear flip armrest with PU armrest pad.
5. Adjustable leg height.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
20KG/22KG	100KG	105CM	59CM	90CM	41CM	40CM	7 inches	10 inches	/	48CM	/	61*36*80CM

YSE143

Ergonomic seat



FEATURES:

Charge time: 4-6H
 Driving Range: 5-10KM
 Max Speed: 6KM/h
 Climbing ability: 9°
 Motor: 250W*2Pcs brushless motor
 Battery: 24V 6Ah/ 10Ah/ 20Ah lithium battery
 Parking Brake: electromagnetic brake
 Cushion: Black

PARAMETERS:

1. Spraying sand black + spraying orange color
2. Comfortable ergonomic seat plate with premium set of removable leather cushions
3. 8 inch front wheel with PU+shock absorber
4. Rear flip armrest with PU armrest pads
5. Leg height adjustable

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
18.5 KG/22 KG	100KG	115CM	59CM	87M	41CM	40CM	8 inches	10 inches	/	48CM	/	61*36*80CM

YSE142

Magnesium alloy



FEATURES:

1. It can be carried on a plane
2. Controller charging
3. Lightweight folding design, easy to carry and travel, easy manual folding
4. Switch to manual model
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15 KM
 Max Speed: 6KM/h
 Climbing ability: 9°
 Motor: 150W*2Pcs brushless motor
 Battery: 24V 6Ah/12Ah lithium
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
18.5/ 22.5KG	100KG	92CM	60CM	93M	43CM	44CM	7 inches	10 inches	48CM	51CM	60*28*70CM	62*35*75CM

YSE143A

Ergonomic seat



FEATURES:

Charge time: 4-6H
 Driving Range: 5-10KM
 Max Speed: 6KM/h
 Climbing ability: 9°
 Motor: 250W*2Pcs
 Battery: 24V 6Ah/ 10Ah/ 20Ah lithium battery
 Parking Brake: electromagnetic brake
 Cushion: Black

PARAMETERS:

1. Spraying sand black
2. Ergonomic second generation backrest, seat plate with advanced set of removable leather cushion.
3. 8 inch front wheel with PU+shock absorber.
4. Rear flip armrest is equipped with PU armrest pad.
5. Adjustable leg height
6. Configuration headrest

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
18.5 KG/22 KG	100KG	115CM	59CM	87M	41CM	40CM	8 inches	10 inches	/	48CM	/	61*36*80CM

Carbon Fiber Electric Wheelchair

YSE133

Carbon fiber



FEATURES:

- 1.It can be carried on a plane
- 2.Lightweight folding design, easy to carry and travel
- 3.Carbon fiber frame material
- 4.Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range:10-15 KM
 Max Speed:8KM/h
 Climbing ability: $\geq 6^\circ$
 Motor: 150W*2Pcs brushless motor
 Battery:24V 8.7Ah/ 15Ah
 Braking distance ≤ 1200 mm
 Obstacle crossing height ≤ 25 mm
 Cushion: gray

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	120KG	102CM	70CM	89CM	45CM	40CM	8 inches	10 inches	42CM	/	70*45*63CM	89*42*80CM

YSE140

Carbon fiber



FEATURES:

- 1.It can be carried on a plane
- 2.Lightweight folding design, easy to carry and travel
- 3.Carbon fiber frame material
- 4.Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range:15 KM
 Max Speed: 6KM/h
 Climbing ability: 9°
 Motor: 220W*2Pcs brushless motor
 Battery:24V 6Ah/12Ah lithium
 Cushion:Black

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	93CM	59CM	92CM	43CM	41CM	7 inches	10 inches	49CM	50CM	59*29*73CM	62*35*75CM

YSE139

Carbon fiber



FEATURES:

- 1.It can be carried on a plane
- 2.Lightweight folding design, easy to carry and travel
- 3.Carbon fiber frame material
- 4.Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range:10-15 KM
 Max Speed:4.5KM/h
 Climbing ability: $\geq 3^\circ$
 Motor: 150W*2Pcs
 Battery:24V 6Ah/10Ah
 Braking distance ≤ 1500 mm
 Obstacle crossing height ≤ 25 mm
 Cushion:Black

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	115CM	96CM	86CM	40CM	40CM	6 inches	12 inches	24CM	43CM	74*36*86CM	89*42*80CM

IYASOCARE®