

IYASOCARE®



FOCUS
RESEARCH
DEVELOPMENT
PRODUCTION

📍 Chongqing Iyasocare Medical Co.,Ltd.
📍 Free Trade Tower, Xiyong street. Chongqing, China
✉️ iyaso@iyasocare.com
☎️ 0086-136 1835 0425


🌐 www.iyasocare.com
iyasocare.en.alibaba.com

CHONGQING IYASOCARE MEDICAL CO.,LTD.

PRODUCT CATALOG BROCHURE

COMPANY INTRODUCTION

Chongqing IYASOCARE MEDICAL CO.,LTD. was established in 2006 and is headquartered in the beautiful mountain city of Chongqing in China. It is an integrated industry and trade company specializing in medical device R&D, production, sales and consulting services, with business all over the world. The company's main products include wheelchairs, hospital beds, electric scooters, walkers, commode chair, etc.; the company is committed to providing customers with reasonable, complete, reliable, and high-quality medical procurement programs; since the company was founded, it has served the United States, the United Kingdom, and Australia, France, Saudi Arabia, Egypt, Pakistan, Japan, South Korea and other global customers; and has obtained medical device registration certificate, medical device production license, medical device business license, ISO13485, CE, TUV, MDR, CFS and other qualifications Certificate; is a supplier of continuous strategic cooperation with China Medical Device Association.



AUTHENTICATION CERTIFICATE

Medical Management System Certification
Award
CHONGQING IYASOCARE MEDICAL CO.,LTD.
Room 305-306, Building 3, Chongqing Railway Post Office Logistics and Storage Project,
Xiyang Group 15th Floor, Xiyang Street, Shapingba District, Chongqing

The Medical management system of the above organizations has been reviewed and complied with **ISO 13485:2016** Requirements of applicable clauses of Medical management system standards

This system is valid for the Production and sales of Medical equipment, Wheelchair, Manual wheelchair, Electric wheelchair, Scooter, Electric handcycle, Manual bed, Electric Bed, Medical bed products

Under the condition that the certificate holder's management system continues to meet the management system standards, the certification deadline is three years. And query on www.ccc-iso.com. The validity of this certificate must be confirmed by the center through regular supervision and review.

Issue Date: December 26, 2023 Valid Period: November 25, 2026

Logos: IAF, UKAS, France, Spain European Certified organization Limited, IVE, etc.

DECLARATION OF CONFORMITY
ACCORDING TO (EU) 2017/745 MEDICAL DEVICE REGULATION

Manufacturer
Name: Chongqing Iyascare Medical Co., Ltd
Address: 406, 4/F, Software Park Phase 2, 8 Xiyuan North Street (Free Trade Tower),
North Street (Free Trade Tower),
The Netherlands
SRN: NL-AR-00000247

Product Information
Name: ELECTRIC WHEELCHAIR
Model: see annex
EMDN: Y122127
Basic UDI DI: /
Classification: Class I, According to Rule 13, Annex VII, Regulation (EU) 2017/745

Declaration
We herewith declare that the above-mentioned products meet the requirements of Medical Device Regulation (EU) 2017/745 and the applicable standards above.

Signature: [Signature] Date: 2023.12.26
Position: GM Place: Chongqing/China

DECLARATION OF CONFORMITY
ACCORDING TO (EU) 2017/745 MEDICAL DEVICE REGULATION

Manufacturer
Name: Chongqing Iyascare Medical Co., Ltd
Address: 406, 4/F, Software Park Phase 2, 8 Xiyuan North Street (Free Trade Tower),
North Street (Free Trade Tower),
The Netherlands
SRN: NL-AR-00000247

Product Information
Name: Electric scooter
Model: see annex
EMDN: Y122124
Basic UDI: /
Classification: Class I, According to Rule 13, Annex VII, Regulation (EU) 2017/745

Declaration
We herewith declare that the above-mentioned products meet the requirements of Medical Device Regulation (EU) 2017/745 and the applicable standards above.

Signature: [Signature] Date: 2023.12.28
Position: GM Place: Chongqing/China

Certificate of Conformity

Certificate No.: ECT2021P095388R-1
Applicant: CHONGQING IYASOCARE MEDICAL CO.,LTD.
Room 201-205,115 Longan Avenue, Tuzhu Town, Shapingba District, Chongqing Province, China
Product Name: Wheelchair, Manual wheelchair, Electric wheelchair, Scooter, Electric handcycle, Commode chair, Walking aids, Transfer wheelchair

Trade Name: N/A
Model: (see the following annex)
Classification: Class I (according to the Manufacturer's declaration)
Review goal: Verification of the presence of Technical Documentation compatible with the Medical Devices Regulation (EU) 2017/745

This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the Technical Documentation shared with us by the manufacturer is compatible with the European Standard for Medical Devices Regulation.

The manufacturer is responsible for the CE Marking process and not exempted to carry out all necessary compliance activities. This document has been issued on the basis of the regulation on ECT Voluntary Mark for the certification of products.

Logos: CE, TÜV Rheinland, Made-in-China

Shenzhen ECT Testing Technology Co., Ltd.
3th Floor, EverBuilding Business Center, Building A, No. 26 Gonghe Industrial Road, Shenzhen, Shenzhen City, China
Tel: +86-755-87587803 Fax: +86-755-87587803 <http://www.ectlab.com>

TÜVRheinland Made-in-China

Report Number: MIC-AS129874 Audit Date: 27 Aug. 2022

This report is issued by Focus Technology Co., Ltd (Made-in-China.com) and the supervising inspectorate (TUV Rheinland) to confirm that:

Company Name: Chongqing Iyascare Medical Co., Ltd.
Showeroom: <https://iyascare.com/made-in-china.com>
Address: Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China
Product: Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products

has been on site audited for the Following Scope of Activity

- General Information
- Foreign Trade Capacity
- Product Research & Development Capacity
- Management System and Product Certification
- Production Capacity & Quality Control
- The actual content of the part relates only to its associated enterprise & Working Environment
- Energy Saving and Emission Reduction
- Photo

General Comments:
Chongqing Iyascare Medical Co., Ltd. was founded in 2010, as a trading company with 10 employees, the company is located at Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China. The company mainly supply Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products. The company has successful foreign trading experience in North America, South America etc.

Signature: [Signature] Date: 16 Feb 2024
On behalf of TÜVRheinland

CEIG Ministerie van Volksgezondheid, Welzijn en Sport

Report Number: MIC-AS129874 Audit Date: 27 Aug. 2022

This report is issued by Focus Technology Co., Ltd (Made-in-China.com) and the supervising inspectorate (TUV Rheinland) to confirm that:

Company Name: Chongqing Iyascare Medical Co., Ltd.
Showeroom: <https://iyascare.com/made-in-china.com>
Address: Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China
Product: Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products

has been on site audited for the Following Scope of Activity

- General Information
- Foreign Trade Capacity
- Product Research & Development Capacity
- Management System and Product Certification
- Production Capacity & Quality Control
- The actual content of the part relates only to its associated enterprise & Working Environment
- Energy Saving and Emission Reduction
- Photo

General Comments:
Chongqing Iyascare Medical Co., Ltd. was founded in 2010, as a trading company with 10 employees, the company is located at Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China. The company mainly supply Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products. The company has successful foreign trading experience in North America, South America etc.

Signature: [Signature] Date: 16 Feb 2024
On behalf of CEIG

Foldable Lightweight Wheelchair

YSM2002



FEATURES:

1. Aluminum frame easy to folding
2. Foldable backrest
3. Rear hand double brakes
4. Flip-up footrest. Backrest with double insurance
5. Leather seat cushion
6. Hydrographics transfer printing carbon fiber surface treatment

Aluminum

YSM2003



FEATURES:

1. Steel frame easy to folding
2. Foldable backrest
3. Rear hand double brakes
4. Flip-up footrest. Backrest with double insurance

Carbon Steel Frame

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
8.9/10KG	100KG	88CM	53CM	89CM	40CM	38CM	6 inches	8 inches	43CM	48CM	72*53*32CM	55*33*73CM

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
11/12KG	100KG	88CM	53CM	89CM	40CM	38CM	6 inches	12 inches	43CM	48CM	72*53*32CM	55*33*73CM

YSM2026



FEATURES:

1. Aluminum alloy frame, sandblack + carbon transfer treatment, rear dual large shock absorbers
2. Ergonomic seat back with removable waterproof cushion.
3. front wheel shock absorber
4. Linkage rear flip armrest with PU armrest pads
5. Self-locking brake
6. Front and rear wheels are equipped with honeycomb non-inflatable shock-absorbing wheels.
7. One button unlocking
8. Adjustable footrest

Aluminum

Ergonomic seat

YSM2027



FEATURES:

1. Aluminum alloy frame, sand black + carbon fiber transfer treatment.
2. Ergonomics II backrest with removable waterproof cushion on the seat back.
3. Front wheel shock absorber
4. Linkage rear flip armrest with PU armrest pads
5. Configure self-locking brake
6. front and rear wheels are equipped with honeycomb non-inflatable shock-absorbing wheels.
7. One-touch unlocking
8. Adjustable footrest
9. Hydraulic reclining adjustment

Aluminum

Ergonomic seat

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
10/12kg	100KG	88CM	55CM	92CM	41CM	40CM	6 inches	12 inches	/	50CM	45*50*73 CM	55*34*82 CM

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
12.5/15kg	100KG	96CM	55CM	98CM	41CM	40CM	11 inches	12 inches	/	50CM	45*50*73 CM	55*34*82 CM

Foldable Lightweight Wheelchair

YSM2028



FEATURES:

1. Aluminum alloy lightweight frame
2. Ergonomic seat backrest, wrapped in mesh black elemental fabric.
3. Aluminum alloy frame with plastic footrests that can be turned up and adjusted.
4. 6-inch wheels, shock absorbing fork. Rear 12-inch wheels
5. With rear hand brake
6. Wheelchair can be folded, armrests can be upturned and fixed, aluminum alloy footrest can be upturned.
7. Telescopic drawbar, easy to pull after folding.
8. Rear push handle is designed.

Aluminum

Ergonomic seat

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	102CM	52.5CM	88CM	42CM	45CM	6 inches	12 inches	43CM	50CM	45*50*73 CM	57*41*85 CM

YSM2029



FEATURES:

1. Aluminum alloy lightweight frame
2. Ergonomic seat backrest, wrapped in mesh black elemental fabric.
3. Aluminum alloy frame with plastic footrests that can be turned up and adjusted.
4. 6-inch honeycomb tires on front wheel, shock absorbing fork. Rear 12-inch honeycomb tires
5. With rear hand brake
6. Wheelchair can be folded, armrests can be upturned and fixed, aluminum alloy footrest can be upturned.
7. Rear push handle is designed and detachable.

Aluminum

Ergonomic seat

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N	100KG	102CM	55CM	88CM	42CM	45CM	6 inches	12 inches	/	46CM	40*55*70 CM	57*39*78 CM

YSM2030



FEATURES:

1. Aluminum alloy frame
2. Light weight
3. Double brakes
4. Flip-up footrest. Backrest with double insurance
5. Flip-up armrest

Aluminum

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N	75KG	74CM	52CM	88CM	34CM	32CM	6 inches	8 inches	41CM	49CM	72*53*32CM	71*30*52CM

IYASOCARE®

Manual Wheelchair - Standard Wheelchair

YSM2001

TYPE: A



FEATURES:

1. Foldable backrest
2. Steel handrail with movable leg folding handle
3. PU handrail, high-grade camouflage 3D breathable seat cushion
4. High quality nylon pedal, protective plate color aluminum plate

Carbon Steel Frame

N/G: 12/14KG

TYPE: B



FEATURES:

1. Aluminum alloy frame light weight
2. Movable handrails, footrest
3. PU armrest, high-grade camouflage 3D breathable, seat cushion, high quality nylon pedal, color aluminum plate
4. 20" aluminum alloy PU tire, equipped with blade light.

Aluminum Frame

N/G: 11/12.8KG

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
/	100KG	95CM	60CM	85CM	44CM	37CM	6 inches	20 inches	65CM	45CM	80*24*70CM	83*26*77CM

YSM2004

TYPE: A



FEATURES:

1. Carbon steel frame, strong and durable
2. Movable handrails, footrest
3. PU armrest, high-grade camouflage 3D breathable, seat cushion, high quality nylon pedal, color aluminum plate

Carbon Steel

N/G: 13/16KG

TYPE: B



FEATURES:

1. Aluminum alloy frame light weight
2. Movable handrails, footrest
3. PU armrest, high-grade camouflage 3D breathable, seat cushion, high quality nylon pedal, color aluminum plate
4. 22" aluminum alloy wheels

Aluminum Frame

N/G: 12/14.8KG

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
/	100KG	104CM	61CM	88CM	46CM	41CM	6 inches	22 inches	/	/	81*26*86CM	82*27*88CM

Manual Wheelchair - Standard Wheelchair

YSM2016

TYPE: A



FEATURES:

- 1.Plating
- 2.Optional pedals: iron, aluminum, plastic
- 3.Armrest optional: blow molding, leather
- 4.Cushion color optional, leather/checkered/black
- 5.Optional tires: solid tires/pneumatic tires

Electroplated

Bedpan

NW: 18.5KG

TYPE: B



FEATURES:

1. Four brakes with toilet
 2. 1.2mm tube thickness
 3. Optional pedals: iron, aluminum, plastic,
 4. Armrest optional: blow molding, leather,
 5. Cushion optional, leather/checkered/black
- Tires: solid tires

Bedpan

NW: 16KG

TYPE: C



FEATURES:

1. Spraying, aluminum alloy rim
2. One potty, mesh breathable cushion
3. Elevating legs, suitable for orthopedics

Orthopedic leg rests

NW: 18.5KG

TYPE: D



FEATURES:

- 1.Comfortable breathable cushion
- 2.Portable handbag, foldable and light for easy travel
- 3.Big rear wheel manual self drive
- 4.Carbon steel raw materials reinforcement design

Integrated wheel

NW: 18.5KG

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	Folding width	DIM
18.5KG	100KG	103CM	68CM	90CM	46CM	45CM	7 inches	24 inches	45CM	52CM	26CM	94*26*93CM

Manual Wheelchair - Standard Wheelchair

YSM2005



FEATURES:

1. Carbon steel frame with gray painted finish
2. reversible armrests, PVC leather armrest pads
3. Detachable footrests, height adjustable.
4. Seat back cushion made of PVC leather, backrest can be folded.
5. 8-inch wider front casters, 24-inch aluminum rear wheels, wider mountain tires
6. Elbow joint brakes, lower than the seat surface after braking.

Offroad

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
20KG/22KG	140KG	100CM	72CM	94CM	56CM	41CM	8 inches	24 inches	40CM	/	100*26*94 CM	81*28*88 CM

YSM2009



FEATURES:

1. PU integral wheel
2. Comfortable honeycomb, breathable cushion
3. Portable handbag, foldable and light for easy travel
4. Four brake designs are smoothly advanced
5. Carbon steel raw materials reinforcement design

Steel

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
12/14KG	100KG	85CM	62CM	82CM	43CM	39CM	6 inches	18 inches	40CM	46CM	85*66*30CM	80*25*70CM

YSM2012



FEATURES:

1. Aluminum alloy two-tone
Frame: Aluminum alloy two-color frame, with frosted sprayed surface.
2. Oxford nylon fabric with mesh cushion, seat back cushion
Seat back pad: Oxford nylon fabric with mesh cushion, seat back pad is made of 600d canvas.
3. Slanted armrests, PU integrated armrest pads.
4. Front wheel rubber tire, aluminum fork.
Rear wheel 20 inches non-inflatable rubber tire, nylon one-piece big wheel.
Large wheels.
5. With rear hand brake function and elbow knuckle brake.
6. Wheelchair back can be folded, armrests can be upturned and fixed, pedals can be upturned.
7. Plastic armrests ring for manual driving of the wheelchair.

Aluminum

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
13/15KG	100KG	91CM	60CM	100CM	45CM	36CM	6 inches	20 inches	/	43CM	91*28*100CM	78*28*70CM

YSM2013



FEATURES:

1. Seat back cushion is made of flame retardant standard oxford nylon fabric mesh cushion, seat back cushion sandwich is 600d canvas.
2. PU one-piece armrest pad.
3. Front wheel 6 inches (15cm)
Rubber tires, high strength steel fork. Rear wheel 20 inches (50cm) solid rubber tire steel wire big wheel.
4. Brake: with rear hand brake function and elbow joint type brake, below the seat surface after braking.
5. Wheelchair back can be folded, armrests can be upturned and fixed, footrests can be upturned. Plastic armrest rings.

Steel

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
12/14KG	100KG	80CM	58CM	90CM	43CM	37CM	6 inches	16 inches	40CM	48CM	78*24*71CM	80*26*73CM

Manual Wheelchair - Standard Wheelchair

YSM2017

TYPE: **A**



FEATURES:

- 1.Fixed armrest and footrest
- 2.Comfortable leather seat cushion
- 3.Foldable and light for easy travel
- 4.Big rear wheel manual self drive
- 5.Carbon steel chromed frame

Plating 809
N/G: 19/21KG

TYPE: **B**



FEATURES:

- 1.Fixed armrest and footrest
- 2.Comfortable leather seat cushion
- 3.Foldable and light for easy travel
- 4.Big rear wheel manual self drive
- 5.Carbon steel chromed frame
- 6.Detachable bedpan
- 7.Safe rear handle brake

Plating 809
Bedpan
N/G: 18/20KG

TYPE: **C**



FEATURES:

1. Net backrest sticky pad
2. Black Oxford cloth cushion
3. Gray PU tire plastic bracelet
4. Plastic pedals
- 5.8 inch narrow front wheel
6. Plastic guard
7. Soft armrests
8. Single brake

Spraying 809
PU rear tire
N/G: 14/16KG

MAIN PARAMETERS												
N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
/	100KG	103CM	68CM	85CM	42CM	45CM	8 inches	24 inches	42CM	51CM	103*27*85CM	94*23*89CM

YSM2032



FEATURES:

1. Aluminum alloy frame, lightweight and fast folding.
2. Quick release rear wheel design, up and down in just one second.
3. Latex cushion, comfortable and worm-proof.
4. Adjustable armrests and removable footrests.
5. Three height adjustments of the rear axle.

Quick release design
Height adjust

MAIN PARAMETERS												
N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
15/17kg	100KG	105CM	70CM	90CM	45CM	42CM	6 inches	12 inches	/	38CM	105*33*73 CM	110*35*76 CM

IYASOCARE®

Manual Wheelchair - Standard Wheelchair

YSM2033



FEATURES:

1. Comfortable honeycomb, breathable cushion
2. Portable handbag, foldable & light for easy travel
3. The widened design makes the riding space bigger and more comfortable
4. Carbon steel raw materials reinforcement design

Wider seat

Spraying

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	Folding size	DIM
19/21KG	100KG	94CM	75CM	110CM	56CM	/	8 inches	24 inches	/	/	79*37*94CM	100*26*90CM

YSM2034



FEATURES:

1. The frame is made of steel tube with spraying treatment, round tube frame design, tube diameter 22mm, wall thickness 1.0mm, the frame is beautiful and solid.
2. The handrail can be upturned to be movable and become a triangular structure, which is convenient for users to get on and off the car, and the handrail pad is made of leather sponge material, which is sturdy and durable.
3. Configuration of plastic footrest, unarmed can be turned back and folded to reduce space for easy storage, footrest and adjustable height, suitable for people of different heights, high strength is not easy to break.

Carbon Steel Frame

Spraying

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding width	DIM
19/20KG	100KG	100CM	68CM	85CM	45CM	45CM	7 inches	24 inches	45CM	49CM	26.5CM	93*25*93CM

YSM2046



FEATURES:

1. Spraying the wheelchair
2. 1 pair of detachable front hanging legs
3. Solid PU wheel
4. The armrest can be folded
5. Full mesh breathable cushion
6. Color optional: turquoise, dark blue, black, red

Carbon Steel Frame

MAIN PARAMETERS

N	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
12KG	200KG	85CM	56CM	85CM	45CM	38CM	6 inches	12 inches	41CM	45CM	/	94*28*89CM

YSM2052



FEATURES:

1. Aluminum alloy painted camouflage one-wheel wheelchair
2. Movable folding handle, 20-inch integrated wheel tire
3. Color: camouflage color

Aluminum Frame

Integrated wheel

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
/	/	95CM	60CM	86CM	44CM	37CM	12CM	30CM	/	45CM	81*60*70CM	70*79*26CM

Manual Wheelchair -Reclining Wheelchair

YSM2021



FEATURES:

1. The frame is made of steel tube with spraying treatment, and it is designed with round tube frame, the diameter of the tube is 22mm, the frame is beautiful and solid.
2. Handrail can not be lifted
3. Configuration of plastic footrest
4. Hanging back is not removable.
5. Leather cushion
6. Rear wheel 24 inch black solid high tyre, front wheel 7 inch solid tyre
7. Mechanical reclining adjustment

Mechanically reclinable

Bedpan

YSM2023



FEATURES:

1. Portable commode
2. Adjustable armrest and footrest
3. Thick and washable bee net
4. Hydraulic control backrest angle
5. Free rubber solid wheel with manual ring
6. Chromed frame, adjustable wheelchair length and height

Hydraulic reclinable

Bedpan

Electroplated

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	Folding width	DIM
N	100KG	104-170CM	68CM	90-140CM	46CM	45CM	7 inches	24 inches	45CM	47CM	26CM	100*26*90CM

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	92-165CM	68CM	90-120CM	45CM	41CM	7 inches	24 inches	45CM	51CM	92*26*92CM	80*30*70CM

YSM2022

TYPE: A



FEATURES:

1. Portable commode
2. Adjustable armrest and footrest
3. Thick and washable bee net
4. Hydraulic control backrest angle
5. Free rubber solid wheel with manual ring
6. Chromed frame, adjustable wheelchair length and height

Hydraulic reclinable

Bedpan

N/G: 20/24KG

rear wheel 24 inches

TYPE: B



FEATURES:

1. Portable commode
2. Adjustable armrest and footrest
3. Thick and washable bee net
4. Hydraulic control backrest angle
5. Free rubber solid wheel with manual ring
6. Chromed frame, adjustable wheelchair length and height

Hydraulic reclinable

Electroplated

N/G: 31/33KG

rear wheel 24 inches

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
/	100KG	92-165CM	68CM	90-120CM	45CM	41CM	7 inches	/	45CM	51CM	92*26*92CM	80*30*70CM

Manual Wheelchairs - Featured Wheelchairs

YSM2025



FEATURES:

1. PU wheels and hand rims.
2. Adjustable belt handhold ring drive.
3. The user's hands do not come into direct contact with the rear wheel, avoiding the risk of 2 cross-infections.
4. More ergonomic, labor-saving and comfortable, to avoid long-term use caused by frozen shoulder and other health hazards.

Ergonomics

Prevention of cross-infection

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	104CM	68CM	93CM	46CM	40CM	8 inches	13 inches	/	45CM	104*37*76 CM	77*40*78 CM

YSM2059



FEATURES:

1. Spray frame, high back, adjustable and detachable headrest
2. Adjustable front and rear seats, adjustable backrest Angle, adjustable armrest height, removable lifting foot support
3. TPR front wheel, plastic wheel frame inflation quickly remove the rear wheel
4. Rear handle with brake

Car seat

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N	150KG	88CM	64CM	121-133CM	46CM	/	8 inches	24 inches	/	58CM	/	/

YSB01



FEATURES:

1. Foldable frame, front and rear wheels quick detachable installation design
2. Adjustable footrest and armrest design, easy to switch the rest position and get on and off the wheelchair.
3. Full steel frame, sturdy and durable, can carry 150kg.
4. Seat and back cushions are made of glass fibre, breathable and waterproof.
5. The front and rear wheels are balloon wheels, which makes it easy to drive on the beach or wade in the water.

Beach

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	136CM	102CM	94CM	46CM	50CM	12* 7 inches	16* 9 inches	32CM	56CM	/	89*42*84CM

YSB02



FEATURES:

1. foldable frame design, front and rear wheels quick release design, can be placed in the trunk of the car.
2. 3-level adjustable leg rest, 3-level adjustable backrest, maximum tilt 150°.
3. Stainless steel frame, waterproof and durable, with a weight capacity of 100 kg.
4. the cushion is made of breathable and waterproof oxford fabric.
5. Front and rear wheels are balloon wheels, easy to drive on the beach or wade in the water.
6. Optional floatation kit for floating in the water.

Beach

Floating

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	120CM	110CM	122CM	50CM	45CM	10 inches	12 inches	/	/	/	82.5*45.5*84.5CM

Manual Wheelchairs - Accessories

Wooden tablecloth		UZ Cushion		16-inch Integrated wheel		20-inch low profile tires	
Plastic dining table		Leather headrest		20-inch Integrated wheel		22-inch low profile tires	
Plastic Hard Handrail		Oxford cloth headrest		24-inch Integrated wheel		24-inch low profile tires	
Leather Foam Armrest		A brake		24-inch nylon Integrated wheel			
Plastic Side Guards		Electroplated hand rim ordinary tire wire wheel		20-inch plastic hand-rimmed wire wheel			
Iron side guards		Normal low tire		22-inch plastic hand-rimmed wire wheel			
Plastic foot pedal		Plain Black Gray Tall Tire		24-inch plastic hand-rimmed wire wheel			
Iron footrest		20-inch pu foam tire		20-inch Electroplated Hand Rim Wire Wheel			
Aluminum foot pedal		24-inch pu foam tires		22-inch Electroplated Hand Rim Wire Wheel			
Fork		24-inch pneumatic tires		24-inch Electroplated Hand Rim Wire Wheel			
Square Basin		6-inch rubber black/ pu gray front wheel		16-inch tires (tall tires)			
Chamber pot		7-inch rubber black/ pu gray front wheel		20-inch tall tires			
Leather leg rests		8-inch rubber black/ pu gray front wheel		22-inch tall tires			
Double Mesh Cushion		12-inch rear wheel		24-inch tall tires			

Sport Wheelchair

YSM2101



FEATURES:

- 1.Fixed baffle
- 2.Folding frame
- 3.Quick remove rear wheel
- 4.Oxford fabric seat and back
- 5.Fixed foot pedal,length adjustable
- 6.Light magnesium alloy frame, surface spray

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
14.5KG/19KG	100KG	80CM	64CM	87CM	38CM	38CM	4 inches	24 inches	34CM	43CM	/	70*63*60CM

YSM2103



FEATURES:

- 1.Frame: the combination of welding of high strength aluminum alloy profile forming, fixed structure, safety performance is good, the surface color spray paint processing, beautiful and durable
- 2.Front wheel: 18 inches PU special racing wheel, adopt the whole aluminum alloy brake, safe and reliable
- 3.Rear wheel: 28 inches inflatable wheels, configuration is 36 knife type high strength aluminum circle, inflatable tires, with light weight, high strength aluminum wheel (hand when using direct drive wheelchair)
- 4.Seat cushion is made of 420D oxford cloth stitching forming, soft texture, high tensile strength, the soft cushion can according to the different requirements of different figure athletes on the wheelchair seat height and the thickness of the configuration to accommodate

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
11.8KG/14KG	100 KG	164CM	63CM	68CM	30-34CM	36CM	18 inches	28 inches	/	45CM	/	165*56*70CM

YSM2102



FEATURES:

- 1.Easy folding, after folding can be put in the trunk
- 2.Quick disassemble rear wheel: Press the quick disassemble button and the rear wheel can be disassembled easily
3. Aviation aluminum alloy frame
4. Anti-overturn wheels are provided on both sides of the back of the wheelchair

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
11.6KG/14KG	100KG	85CM	56CM	77CM	36CM	48CM	4 inches	24 inches	/	48CM	85*34*77CM	87*36*79CM

YSM2103C



FEATURES:

- 1.Frame: the combination of welding of high strength aluminum alloy profile forming, fixed structure, safety performance is good, the surface color spray paint processing, beautiful and durable
- 2.Front wheel: 18 inches PU special racing wheel, adopt the whole aluminum alloy brake, safe and reliable
- 3.Rear wheel: 28 inches inflatable wheels, configuration is 36 knife type high strength aluminum circle, inflatable tires, with light weight, high strength aluminum wheel (hand when using direct drive wheelchair)
- 4.Seat cushion is made of 420D oxford cloth stitching forming, soft texture, high tensile strength, the soft cushion can according to the different requirements of different figure athletes on the wheelchair seat height and the thickness of the configuration to accommodate

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
12KG/15KG	100KG	164CM	63CM	68CM	28/30/32 CM	36CM	18 inches	26 inches	/	45CM	/	165*65*70CM

Sport Wheelchair

YSM2104



FEATURES:

1. Easy folding, after folding can be put in the trunk
2. Quick disassemble rear wheel: Press the quick disassemble button and the rear wheel can be disassembled easily
3. Aviation aluminum alloy frame, designed for dance moves
4. Anti-overturn wheels are provided on both sides of the back of the wheelchair

Dance

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	85CM	56CM	75CM	36CM	/	2.4 inches	22 inches	/	50CM	/	85*58*80CM

YSM2106



FEATURES:

1. Easy folding, after folding can be put in the trunk
2. Quick disassemble rear wheel: Press the quick disassemble button and the rear wheel can be disassembled easily
3. Aviation aluminum alloy frame, designed for tennis
4. Anti-overturn wheels are provided on both sides of the back of the wheelchair

Tennis

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	85CM	85CM	68CM	30/32/34/36/38/40CM	/	3.2 inches	26 inches	/	42CM	/	103*51*67CM

YSM2105



FEATURES:

1. Easy folding, after folding can be put in the trunk
2. Quick disassemble rear wheel: Press the quick disassemble button and the rear wheel can be disassembled easily
3. Aviation aluminum alloy frame, designed for basketball guard
4. Anti-overturn wheels are provided on both sides of the back of the wheelchair

Basketball

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	85CM	85CM	64CM	30/32/34/36/38/40CM	/	4 inches	26 inches	/	50CM	/	82*80*86CM

YSM2107



FEATURES:

1. Aluminum wheelchair frame
2. Quick release of rear wheel, convenient storage of wheelchair
3. Wear a seat belt
4. Sturdy plastic footboards

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	75CM	67CM	77CM	36CM	/	4 inches	24 inches	/	48CM	/	76*62*50CM

Sport Wheelchair

YSM2108



FEATURES:

1. Easy folding, after folding can be put in the trunk
2. Quick disassemble rear wheel: Press the quick disassemble button and the rear wheel can be disassembled easily
3. Aviation aluminum alloy frame
4. Anti-overturn wheels are provided on both sides of the back of the wheelchair

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	81.5CM	57-63CM	76.5CM	36/38/40/42CM	40CM	4 inches	24 inches	/	46CM	81.5*35*76.5CM	69*38*85CM

YSM2111



FEATURES:

1. Lightweight design, fashion modeling
2. Carbon fiber frame
3. Anti-overturn wheels are provided on both sides of the back of the wheelchair

Carbon fiber

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	84CM	63-69CM	77CM	36/38/40/42CM	40CM	4 inches	24 inches	/	46CM	/	90*65*78CM

YSM2109



FEATURES:

1. Aluminum foldable frame
2. Fixed armrests, foldable and height-adjustable footrests
3. wider solid front wheels, wider off-road quick release light-weight pneumatic rear wheels
4. Comfortable foam seat and back cushion

All terrain

Offroad

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	84CM	60CM	77CM	36/38/40/42CM	37CM	4 inches	24 inches	/	48CM	87*33*76 CM	87*36*79 CM

YSM2112



FEATURES:

1. Lightweight design, fashion modeling
2. Aluminum frame
3. Aluminum alloy manual parking brake
4. Cushion fabric selection of high-quality mesh fabric lined with high-quality high resilience sponge

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
16/18KG	100KG	80CM	54CM	76CM	36/38/40/42CM	37CM	5 inches	24 inches	30CM	52CM	/	68*64*58CM

Cerebral Palsy Wheelchair

YSM2201



FEATURES:

1. One-piece Oxford cloth backrest cushion with butterfly seat belt and double brake;
2. ABS guard plate. Feet can be lifted up and removable. Easy to carry;
3. ABS foot pedal, PU handrail, upper and lower three gears adjustable.
4. 135° adjustable back, 90° backrest, adjustable Angle.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
22KG/24KG	100KG	112CM	46CM	100CM	38CM	39CM	6 inches	16 inches	/	46CM	80*25*90CM	83*30*94CM

YSM2205



FEATURES:

1. Butterfly seat belt and double brake;
2. Feet can be lifted up and removable. Easy to carry;
3. Aluminum alloy front fork, ABS foot pedal, PU handrail, upper and lower three gears adjustable.
4. 150° adjustable back, 90° backrest, adjustable Angle.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
18/22KG	75KG	92CM	42CM	103CM	30CM	30CM	6 inches	12 inches	45CM	49CM	75*37*90CM	80*41*69CM

YSM2201L



FEATURES:

1. One-piece Oxford cloth backrest cushion with butterfly seat belt and double brake;
2. ABS guard plate. Feet can be lifted up and removable. Easy to carry;
3. ABS foot pedal, PU handrail, upper and lower three gears adjustable.
4. 135° adjustable back, 90° backrest, adjustable Angle.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
18KG/20KG	100KG	114CM	46CM	100CM	38CM	35CM	6 inches	16 inches	45CM	50CM	84*32*95CM	83*30*94CM

YSM2206



FEATURES:

1. Aluminum alloy spray frame;
2. 3D mesh surface cushion, PVC wheels;
3. 150° adjustable back, 90° backrest, adjustable Angle.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
25KG/27KG	75KG	95CM	46CM	112CM	33CM	35CM	6 inches	12 inches	45CM	60-72CM	75*37*90CM	85*51*107CM

Children Wheelchair

YSE402

Electric



FEATURES:

1. Aluminum frame lightweight folding design, easy to carry and travel
2. Smart controller, easy to study
3. Detachable footrest
4. Double back handle
5. Flip-up armrest

PARAMETERS:

Climb Angle: 8 degrees
 Motor: 200w*2 brushless motor
 Battery: 24V 10Ah lithium
 Charge time: 4-6 H
 Maximum speed: 6 km/h
 Maximum range: 15 km

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	100CM	52CM	102CM	39CM	39CM	7 inches	12 inches	28CM	58CM	74*31*69 CM	79*38*75 CM

YSM2202



FEATURES:

1. Aluminum alloy frame, tube diameter 25MM, wall thickness 1.5MM, baking paint.
2. the footrest height can be adjusted and can be rotated to facilitate easy to get on and off the bike, the whole bike can be folded
3. The seat is high in the front and low in the back, and equipped with a strap-type three-point safety belt to prevent the child from slipping down.
4. 16 inch PU rear wheels, aluminum alloy parking brake.
5. With sunshade, with armrests handrail

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	115CM	57CM	107CM	44CM	35-45CM	8 inches	16 inches	62CM	55-45CM	120*38*35CM	120*38*35CM

YSE403

Electric



FEATURES:

1. Aluminum frame lightweight folding design, easy to carry and travel
2. Smart controller, easy to study
3. Detachable footrest
4. Double back handle
5. Flip-up armrest

PARAMETERS:

Climb Angle: 8 degrees
 Motor: 200w*2 worm gear motors
 Battery: 24V 23Ah lithium
 Charge time: 4-6 H
 Maximum speed: 6 km/h
 Maximum range: 20-25 km

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100kg	96CM	56CM	92CM	38CM	40CM	7 inches	12 inches	/	52CM	74*31*69CM	86*58*75CM

YSM2203



FEATURES:

1. The surface is treated with baking paint, which has the characteristics of no fading and strong anti-rust ability.
2. Foldable structure with high safety performance
3. Oxford cloth fabric sewing molding, soft texture, breathable and non-slip
4. backrest armrest can be adjusted, footrest can be detached.

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	100CM	51CM	89CM	35CM	36CM	6 inches	22 inches	41CM	49CM	/	78*41*92CM

Children Wheelchair

YSM2204



FEATURES:

1. An adaptable aluminium manual wheelchair designed for children with disabilities.
2. The seat width and seat depth can be adjusted. With the growth of children, the seat width can be adjusted to 25cm, 30cm and 35cm
3. Seat depth can be adjusted to 25cm, 30cm, 35cm three kinds of seat depth.

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	75KG	92CM	60CM	92CM	25-35CM	25-35CM	4 inches	20 inches	/	36-41CM	92*36*87CM	78*41*92CM

YSM2209



FEATURES:

1. Comfortable and breathable cushion
2. Aluminium alloy frame
3. Solid rubber wheels with manual ring
4. Adjustable handrest
5. Removable footrest

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	75KG	95CM	46CM	90CM	35CM	36CM	6/ inches	22/ inches	38CM	48CM	83*28*90CM	85*29*92CM

YSM2208



FEATURES:

1. Comfortable and breathable cushion
2. Aluminium alloy frame
3. Solid rubber wheels with manual ring
4. Rear brake handle
5. Removable footrest

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	75KG	75CM	48CM	94CM	30CM	30CM	4 inches	20 inches	38CM	46CM	65*30*69CM	71*33*69CM

YSM-PONY



FEATURES:

1. the frame aluminum alloy material, the surface is painted pearl white, the whole vehicle can be folded.
2. backrest 3 gear can be adjusted back 90-110 degrees, 5 inch universal front wheel, 5 inch directional rear wheel, treadle brake.
3. 3 gear sitting depth can be adjusted.
4. foot support height is adjustable.
5. the head pillow height is adjustable.
6. hips with 3-point seat belt.
7. maximum load 75KG.

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	75KG	96CM	40CM	94CM	33CM	25CM	/ inches	/ inches	/	48CM	103*40*34CM	110*47*42CM

Shower And Commode Wheelchair

YSM2007



FEATURES:

- 1, Three-in-one wheelchair, commode chair, shower chair, manual wheelchair
- 2, Detachable commode for easy cleaning
- 3, Detachable footrest, easy to fold
- 4, Adjustable leather armrests
- 5, Four small wheels with independent brakes

YSM2057



FEATURES:

1. Aluminum alloy frame with folding handle
2. Fabric material: luxurious sofa cushion
3. Wheel Material: Shock Absorbing PU
4. Waterproofing

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
13/15KG	100KG	100CM	54CM	93CM	41CM	41CM	/	/	/	47CM	81*26*86CM	67*60*66CM

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
13/15KG	100KG	100CM	54CM	93CM	41CM	41CM	/	/	/	47CM	81*26*86CM	67*60*66CM

YSM2058



FEATURES:

1. Aluminum alloy frame with folding handle
2. Fabric material: luxurious sofa cushion
3. Wheel Material: Shock Absorbing PU
4. Waterproofing

MAIN PARAMETERS

NW	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
13KG	100 KG	91CM	57CM	90CM	/	/	6 inches	20 inches	40CM	/	91*25*70CM	71*29*70CM

IYASOCARE®