

IYASOCARE®



FOCUS
RESEARCH
DEVELOPMENT
PRODUCTION

📍 Chongqing Iyasocare Medical Co.,Ltd.
📍 Free Trade Tower, Xiyong street. Chongqing, China
✉️ iyaso@iyasocare.com
☎️ 0086-136 1835 0425


🌐 www.iyasocare.com
iyasocare.en.alibaba.com

CHONGQING IYASOCARE MEDICAL CO.,LTD.

PRODUCT CATALOG BROCHURE

COMPANY INTRODUCTION

Chongqing IYASOCARE MEDICAL CO.,LTD. was established in 2006 and is headquartered in the beautiful mountain city of Chongqing in China. It is an integrated industry and trade company specializing in medical device R&D, production, sales and consulting services, with business all over the world. The company's main products include wheelchairs, hospital beds, electric scooters, walkers, commode chair, etc.; the company is committed to providing customers with reasonable, complete, reliable, and high-quality medical procurement programs; since the company was founded, it has served the United States, the United Kingdom, and Australia, France, Saudi Arabia, Egypt, Pakistan, Japan, South Korea and other global customers; and has obtained medical device registration certificate, medical device production license, medical device business license, ISO13485, CE, TUV, MDR, CFS and other qualifications Certificate; is a supplier of continuous strategic cooperation with China Medical Device Association.



AUTHENTICATION CERTIFICATE

Medical Management System Certification
Award
CHONGQING IYASOCARE MEDICAL CO.,LTD.
Room 305-306, Building 8, Chongqing Railway Post Office Logistics and Storage Project,
Xiyang Group 15th Floor, Xiyang Street, Shapingba District, Chongqing

The Medical management system of the above organizations has been reviewed and complied with **ISO 13485:2016** Requirements of applicable clauses of Medical management system standards

This system is valid for the Production and sales of Medical equipment, Wheelchair, Manual wheelchair, Electric wheelchair, Scooter, Electric handcycle, Manual bed, Electric Bed, Medical bed products

Under the condition that the certificate holder's management system continues to meet the management system standards, the certification deadline is three years. And query on www.ccc-iso.com. The validity of this certificate must be confirmed by the center through regular supervision and review.

Issue date: December 26, 2023 Valid period: November 25, 2028

Spain European Certified organization Limited
www.ccc-iso.com

CHONGQING IYASOCARE MEDICAL CO.,LTD.
重庆伊亚索凯医疗器械有限公司

DECLARATION OF CONFORMITY
ACCORDING TO (EU) 2017/745 MEDICAL DEVICE REGULATION

Manufacturer
Name: Chongqing Iyascare Medical Co., Ltd
Address: 406, 4/F, Software Park Phase 2, 8 Xiyuan North Street (Free Trade Tower),
The Netherlands
SRN: NL-AR-00000247

Product Information
Name: ELECTRIC WHEELCHAIR
Model: see annex
EMDN: Y122127
Basic UDI DI: /
Classification: Class I, According to Rule 13, Annex VII, Regulation (EU) 2017/745

Declaration
We herewith declare that the above-mentioned products meet the requirements of Medical Device Regulation (EU) 2017/745 and the applicable standards above.

Signature: [Signature] Date: 2023.12.28
Position: GM Place: Chongqing/China

DECLARATION OF CONFORMITY
ACCORDING TO (EU) 2017/745 MEDICAL DEVICE REGULATION

Manufacturer
Name: Chongqing Iyascare Medical Co., Ltd
Address: 406, 4/F, Software Park Phase 2, 8 Xiyuan North Street (Free Trade Tower),
The Netherlands
SRN: NL-AR-00000247

Product Information
Name: Electric scooter
Model: see annex
EMDN: Y122124
Basic UDI: /
Classification: Class I, According to Rule 13, Annex VII, Regulation (EU) 2017/745

Declaration
We herewith declare that the above-mentioned products meet the requirements of Medical Device Regulation (EU) 2017/745 and the applicable standards above.

Signature: [Signature] Date: 2023.12.28
Position: GM Place: Chongqing/China

Certificate of Conformity

Certificate No.: ECT2021P095388R-1
Applicant: CHONGQING IYASOCARE MEDICAL CO.,LTD.
Room 201-205,115 Longan Avenue, Tuzhu Town, Shapingba District, Chongqing Province, China
Product Name: Wheelchair, Manual wheelchair, Electric wheelchair, Scooter, Electric handcycle, Commode chair, Walking aids, Transfer wheelchair

Trade Name: N/A
Model: (see the following annex)
Classification: Class I (according to the Manufacturer's declaration)
Review goal: Verification of the presence of Technical Documentation compatible with the Medical Devices Regulation (EU) 2017/745

This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the Technical Documentation shared with us by the manufacturer is compatible with the European Standard for Medical Devices Regulation.

The manufacturer is responsible for the CE Marking process and not exempted to carry out all necessary compliance activities. This document has been issued on the basis of the regulation on ECT Voluntary Mark for the certification of products.

Shenzhen ECT Testing Technology Co., Ltd.
3th Floor, EverBuilding Business Center, Building A, No. 26 Gonghe Industrial Road, Shenzhen, Shenzhen City, China.
Tel: +86-755-87587803 Fax: +86-755-87587803 <http://www.ectlab.com>

TÜVRheinland Made-in-China

Report Number: MIC-AS129874 Audit Date: 27 Aug. 2022

This report is issued by Focus Technology Co., Ltd (Made-in-China.com) and the supervising inspectorate (TUV Rheinland) to confirm that:

Company Name: Chongqing Iyascare Medical Co., Ltd.
Showeroom: <https://iyascare.com/made-in-china>
Address: Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China
Product: Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products

has been on site audited for the Following Scope of Activity

- General Information
- Foreign Trade Capacity
- Product Research & Development Capacity
- Management System and Product Certification
- Production Capacity & Quality Control
- The actual content of the part relates only to its associated enterprise & Working Environment
- Energy Saving and Emission Reduction
- Photo

General Comments:
Chongqing Iyascare Medical Co., Ltd. was founded in 2010, is a trading company with 10 employees, the company is located at Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China. The company mainly supply Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products. The company has successful foreign trading experience in North America, South America etc.

Signature: [Signature] Date: 16 Feb 2024
TUV Rheinland

CEIG Ministerie van Volksgezondheid, Welzijn en Sport

Report Number: MIC-AS129874 Audit Date: 27 Aug. 2022

This report is issued by Focus Technology Co., Ltd (Made-in-China.com) and the supervising inspectorate (TUV Rheinland) to confirm that:

Company Name: Chongqing Iyascare Medical Co., Ltd.
Showeroom: <https://iyascare.com/made-in-china>
Address: Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China
Product: Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products

has been on site audited for the Following Scope of Activity

- General Information
- Foreign Trade Capacity
- Product Research & Development Capacity
- Management System and Product Certification
- Production Capacity & Quality Control
- The actual content of the part relates only to its associated enterprise & Working Environment
- Energy Saving and Emission Reduction
- Photo

General Comments:
Chongqing Iyascare Medical Co., Ltd. was founded in 2010, is a trading company with 10 employees, the company is located at Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China. The company mainly supply Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products. The company has successful foreign trading experience in North America, South America etc.

Signature: [Signature] Date: 16 Feb 2024
TUV Rheinland

Standing Wheelchair

YSE301

Full lying



FEATURES:

1. Electromagnetic brake system, carbon steel body,
2. Hydraulic adjustment device can be completely lay down.
3. Front and rear shock absorption design, rocker anti-radio interference controller.
4. Equipped with high-quality knee pads.

PARAMETERS:

Driving Range: 45KM
 Max Speed: 8KM/h
 Battery: 24V 50A(lead acid)
 Motor: 24V 350W*2
 Parking Brake: Intelligent electronic brake + electromagnetic brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	120CM	65CM	120CM	48CM	45CM	14 inches	10 inches	52CM	55CM	/	125*70*70CM

YSE303



FEATURES:

1. Electric walking, electric standing, electric standing and walking
2. Turn on the spot with 360 degrees
3. One-hand joystick controller
4. Universal front wheel
5. Armrests on both sides can be raised

PARAMETERS:

Driving Range: 20-25 KM
 Max Speed: 7KM/h
 Battery: 24V 28A(lead acid)
 Motor: 24V 200W*2 worm gear motor
 Parking: Electronic brake
 Cushion: Black with red

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	120KG	120CM	65CM	110CM	45CM	45CM	8 inches	15 inches	/	55CM	/	112*72*95CM

YSE302

All terrain



FEATURES:

1. Suitable for disabled and elderly people's electric powerful cross-country wheelchair.
2. It can climb stairs of 10 centimeters high.
3. Anti-slip slope design, the wheelchair will not slide backward or roll over on the slope.
4. Comes standard with USB charging port, leather seat, seat belt, turn signal, headlight and horn.
5. Speed and battery life vary with battery and motor power.
6. Three driving modes are available: off-road, normal, and indoor.
7. The seat can be adjusted to 180 degrees of recline.
8. Motorized standing available with remote control to assist standing anytime, anywhere.

PARAMETERS:

Charging time: 8-10h
 Driving range :40-60km
 Top speed :28 km/h
 Climbing capacity :15°
 Battery: 48v 20ah lead acid
 Parking brake: electromagnetic brake
 Cushion:black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	200KG	125CM	64CM	108CM	49CM	46CM	14 inches	14 inches	/	/	/	135*70*90CM

YSE304

Full lying



FEATURES:

1. Electromagnetic brake system, carbon steel body.
2. It can be electrically adjusted to standing mode and full lying mode, manual electric dual-mode system, multi-function remote control operating panel, each function button exists independently, optional toggle switch, operated by wrist.
3. Large-capacity battery, long life use.
4. Equipped with high quality knee pads and double-sided lifting armrests.

PARAMETERS:

Driving Range: 45KM
 Max Speed: 8KM/h
 Battery: 24V 55Ah
 Motor: 24V 450W
 Parking Brake: Intelligent electronic brake + electromagnetic brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	120KG	120CM	65CM	110CM	45CM	45CM	8 inches	15 inches	/	55CM	/	112*72*95CM

Standing Wheelchair

YSE305



FEATURES:

1. Electric walking, electric standing, electric standing and walking
2. Turn on the spot with 360 degrees
3. One-hand joystick controller
4. Universal front wheel
5. Armrests on both sides can be raised

PARAMETERS:

Driving Range: 30-40 KM
 Max Speed: 8KM/h
 Battery: 24V 40A (lead acid)
 Motor: 250W*2 24V
 Parking Brake: Electromagnetic Brake and intelligent electronic brake
 Cushion: Black with red

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	100CM	56CM	85M	45CM	43CM	8 inches	14/22 inches	/	47CM	/	115*65*60CM

YSE308

Full lying



FEATURES:

1. Electromagnetic brake system, carbon steel body.
2. 4 modes: seating/elevating/reclining/standing. Multi-function remote control operation panel, each function button is independent, optional toggle switch, operated by wrist.
3. Large capacity battery with long service life.
4. Equipped with high quality knee pads and double-sided lifting armrests.

PARAMETERS:

Driving Range: 45KM
 Max Speed: 6KM/h
 Gradeability: 12°
 Battery: 24V 55Ah
 Motor: 320w*2 worm gear motor
 Parking Brake: Intelligent electronic brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	120KG	132CM	65CM	133M	43CM	46CM	8 inches	14 inches	54CM	56-69CM	/	118*71*86CM

YSE306

All terrain



FEATURES:

Charging time: 6-8h
 Driving range :25-60km
 Motor: 48V 600W motor
 Top speed :20 km/h
 Climbing capacity :15°
 Battery: 48v 20ah lead acid
 Parking brake: drum brake
 Cushion: black/ red

PARAMETERS:

1. Electric off-road wheelchair suitable for disabled and elderly people.
2. Slope anti-slip design
3. Standard equipped with standing controller
4. 3 speed adjustable, maximum speed 20km/h. 8.
5. Motorized standing by remote control

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	120KG	125CM	65CM	125M	45CM	49CM	14 inches	13 inches	/	53CM	/	125*70*68CM

YSE309



FEATURES:

1. Electromagnetic brake system, carbon steel body,
2. Hydraulic adjustment device can be completely lay down.
3. Front and rear shock absorption design, rocker anti-radio interference controller.
4. Equipped with high-quality knee pads, and double-sided lifting armrests.

PARAMETERS:

Driving Range: 8-15km
 Max Speed: 0-8 km/h
 Battery: 20A lithium battery
 Parking Brake: hand brake + electromagnetic brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	108CM	68CM	130CM	49CM	45CM	12 inches	8 inches	/	62CM	/	/

Standing Wheelchair

YSE310

Rehabilitation exercise



FEATURES:

- 1.Independent standing function, sitting/standing switching
- 2.Auxiliary walking training, correcting the patient's spine
- 3.Seamless aluminum alloy frame, lightweight design
- 4.Electromagnetic brake, safe parking

PARAMETERS:

Driving Range:20km
 Max Speed:0-6 km/h
 Gradeability: 10°
 Battery: 24V 15.4A lithium
 Parking Brake: electromagnetic brake
 Cushion:Black/gray

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	105CM	69CM	109CM	49CM	45CM	7 inches	10 inches	/	/	/	/

YSM4005

Manual standing



FEATURES:

- 1.Manual standing up
- 2.Comfortable cushion
- 3.Free rubber big rear wheel with manual ring
- 4.Steel frame,adjustable wheelchair height

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	106CM	67CM	85CM	42CM	43CM	5 inches	24 inches	/	45CM	/	115*65*64CM

YSE311

Semi-automatic standing



FEATURES:

Power standing wheelchair, armrest, comfortable seat belt, power control standing, manual control of direction and progress.

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100 KG	/	65CM	104CM	46CM	/	8 inches	24 inches	/	58CM	/	110*69*67CM

IYASOCARE®

Electric Lifting Wheelchair

YSE307



FEATURES:

1. One-hand joystick controller, front wheel universal steering
2. The armrests on both sides can be quickly removed, the backrest can be reclined, the pedal can be quickly removed, and the lifting pedal can be replaced
3. Turned on the spot in a 360-degree,
4. It has electric walking and electric lifting functions

PARAMETERS:

Driving Range: 25 km
 Max Speed: 0-8 km/h
 Motor: 24V 250W*2
 Battery: 24V 24A lead acid
 Parking Brake: Intelligent electronic brake + electromagnetic brake
 Cushion: Black and red

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	105CM	59CM	110-165CM	45CM	45CM	8 inches	12 inches	72CM	45-95CM	/	/

YSE314



FEATURES:

1. One-hand joystick controller, front wheel universal steering
2. Luxurious and comfortable car seats for long periods of use
3. It has electric driving and electric lifting functions

PARAMETERS:

Charge time: 8-10H
 Driving Range: 15 KM
 Max Speed: 6KM/h
 Climbing ability: 12°
 Motor: 24V 180W*2
 Battery: 20AH lithium
 Brake type: Electromagnetic brake
 Braking distance ≤ 800mm
 Obstacle crossing height ≤ 40mm
 Cushion: Black

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	125KG	93CM	59CM	140CM	47CM	45CM	8 inches	12.5 inches	/	/	75*59*69CM	77*63*72KG

YSE313



FEATURES:

1. Universal controller
2. Electric strong power offroad wheels wheelchair for disabled and elderly people.
3. Anti-slip design, the wheelchair will not slide backwards on the slope
4. Safe low-speed start design to avoid potential safety hazards caused by excessive acceleration
5. The seat can go up 45 degrees, and the back can lie down 160 degrees with electric adjustment.
6. The foot rest can go up 70 degrees adjustable.
7. Seat rotation design, 360° seat free adjustment.
8. Lift seat and assisted standing UP

PARAMETERS:

Charge time: 8-10H
 Driving Range: 20-40 KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Battery: 24V 10.4Ah/15.6Ah lithium
 Brake type: Electromagnetic brake
 Braking distance ≤ 800mm
 Obstacle crossing height ≤ 80mm
 Cushion: Black

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	92CM	68CM	115CM	48CM	42CM	12 inches	8 inches	23CM	47CM	/	/

YSE315



FEATURES:

1. One-handed joystick controller
2. Adjustable armrests and footrests on both sides.
3. 360 degrees flip in place
4. Remote control lift, seat front and rear adjustable, backrest adjustable

PARAMETERS:

Charge time: 4-6H
 Driving Range: 15-25 KM
 Max Speed: 8KM/h
 Climbing ability: 12°
 Motor: 250W*2
 Battery: 12/ 15/ 20AH
 Brake type: Electromagnetic brake
 Braking distance ≤ 1500mm
 Cushion: Black

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	120KG	102CM	65CM	105CM	42CM	40CM	8 inches	12 inches	/	43CM	/	77*60*85CM