

IYASOCARE®



FOCUS
RESEARCH
DEVELOPMENT
PRODUCTION

📍 Chongqing Iyasocare Medical Co.,Ltd.
📍 Free Trade Tower, Xiyong street. Chongqing, China
✉️ iyaso@iyasocare.com
☎️ 0086-136 1835 0425


🌐 www.iyasocare.com
iyasocare.en.alibaba.com

CHONGQING IYASOCARE MEDICAL CO.,LTD.

PRODUCT CATALOG BROCHURE

COMPANY INTRODUCTION

Chongqing IYASOCARE MEDICAL CO.,LTD. was established in 2006 and is headquartered in the beautiful mountain city of Chongqing in China. It is an integrated industry and trade company specializing in medical device R&D, production, sales and consulting services, with business all over the world. The company's main products include wheelchairs, hospital beds, electric scooters, walkers, commode chair, etc.; the company is committed to providing customers with reasonable, complete, reliable, and high-quality medical procurement programs; since the company was founded, it has served the United States, the United Kingdom, and Australia, France, Saudi Arabia, Egypt, Pakistan, Japan, South Korea and other global customers; and has obtained medical device registration certificate, medical device production license, medical device business license, ISO13485, CE, TUV, MDR, CFS and other qualifications Certificate; is a supplier of continuous strategic cooperation with China Medical Device Association.



AUTHENTICATION CERTIFICATE

Medical Management System Certification
Award
CHONGQING IYASOCARE MEDICAL CO.,LTD.
Room 305-306, Building 3, Chongqing Railway Post Office Logistics and Storage Project, Xiyang Group 1# Standard Plaza, Yuzhou Street, Shapingba District, Chongqing

The Medical management system of the above organizations has been reviewed and complied with **ISO 13485:2016** Requirements of applicable clauses of Medical management system standards

This system is valid for the Production and sales of Medical equipment, Wheelchair, Manual wheelchair, Electric wheelchair, Scooter, Electric handcycle, Manual bed, Electric Bed, Medical bed products

Under the condition that the certificate holder's management system continues to meet the management system standards, the certification deadline is three years. And query on www.ccc-iso.com. The validity of this certificate must be confirmed by the center through regular supervision and review.

Issue date: December 26, 2023 Valid period: November 25, 2026

Spain European Certified organization Limited
www.ccc-iso.com

China Medical Device Association
www.cmdea.com

DECLARATION OF CONFORMITY
ACCORDING TO (EU) 2017/745 MEDICAL DEVICE REGULATION

Manufacturer
Name: Chongqing Iyascare Medical Co., Ltd
Address: 406, 4/F, Software Park Phase 2, 8 Xiyuan North Street (Free Trade Tower), The Netherlands
SRN: NL-AR-00000247

Product Information
Name: ELECTRIC WHEELCHAIR
Model: see annex
EMDN: Y122127
Basic UDI DI: /
Classification: Class I, According to Rule 13, Annex VII, Regulation (EU) 2017/745

Declaration
We herewith declare that the above-mentioned products meet the requirements of Medical Device Regulation (EU) 2017/745 and the applicable standards above.

Signature: [Signature] Date: 2023.12.26
Position: GM Place: Chongqing/China

DECLARATION OF CONFORMITY
ACCORDING TO (EU) 2017/745 MEDICAL DEVICE REGULATION

Manufacturer
Name: Chongqing Iyascare Medical Co., Ltd
Address: 406, 4/F, Software Park Phase 2, 8 Xiyuan North Street (Free Trade Tower), The Netherlands
SRN: NL-AR-00000247

Product Information
Name: Electric scooter
Model: see annex
EMDN: Y122124
Basic UDI: /
Classification: Class I, According to Rule 13, Annex VII, Regulation (EU) 2017/745

Declaration
We herewith declare that the above-mentioned products meet the requirements of Medical Device Regulation (EU) 2017/745 and the applicable standards above.

Signature: [Signature] Date: 2023.12.28
Position: GM Place: Chongqing/China

Certificate of Conformity

Certificate No.: ECT2021P095388R-1
Applicant: CHONGQING IYASOCARE MEDICAL CO.,LTD.
Room 201-205,115 Longan Avenue, Tuzhu Town, Shapingba District, Chongqing Province, China
Product Name: Wheelchair, Manual wheelchair, Electric wheelchair, Scooter, Electric handcycle, Commode chair, Walking aids, Transfer wheelchair

Trade Name: N/A
Model: (see the following annex)
Classification: Class I (according to the Manufacturer's declaration)
Review goal: Verification of the presence of Technical Documentation compatible with the Medical Devices Regulation (EU) 2017/745

This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the Technical Documentation shared with us by the manufacturer is compatible with the European Standard for Medical Devices Regulation.

The manufacturer is responsible for the CE Marking process and not exempted to carry out all necessary compliance activities. This document has been issued on the basis of the regulation on ECT Voluntary Mark for the certification of products.

Shenzhen ECT Testing Technology Co., Ltd.
3th Floor, EverBuilding Business Center, Building A, No. 26 Gonghe Industrial Road, Shenzhen, Shenzhen City, China
Tel: +86-755-87587803 Fax: +86-755-87587803 <http://www.ectlab.com>

TÜVRheinland Made-in-China

Report Number: MIC-AS129874 Audit Date: 27 Aug. 2022

This report is issued by Focus Technology Co., Ltd (Made-in-China.com) and the supervising inspectorate (TUV Rheinland) to confirm that:

Company Name: Chongqing Iyascare Medical Co., Ltd.
Showroom: <https://iyascare.com/made-in-china.com>
Address: Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China
Product: Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products

has been on site audited for the following Scope of Activity

- General Information
- Foreign Trade Capacity
- Product Research & Development Capacity
- Management System and Product Certification
- Production Capacity & Quality Control
- The actual content of the part relates only to its associated enterprise & Working Environment
- Energy Saving and Emission Reduction
- Photo

General Comments:
Chongqing Iyascare Medical Co., Ltd. was founded in 2010, is a trading company with 10 employees, the company is located at Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China. The company mainly supply Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products. The company has successful foreign trading experience in North America, South America etc.

Signature: [Signature] Date: 16 Sep 2021

CEIG Ministerie van Volksgezondheid, Welzijn en Sport

Report Number: MIC-AS129874 Audit Date: 27 Aug. 2022

This report is issued by Focus Technology Co., Ltd (Made-in-China.com) and the supervising inspectorate (TUV Rheinland) to confirm that:

Company Name: Chongqing Iyascare Medical Co., Ltd.
Showroom: <https://iyascare.com/made-in-china.com>
Address: Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China
Product: Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products

has been on site audited for the following Scope of Activity

- General Information
- Foreign Trade Capacity
- Product Research & Development Capacity
- Management System and Product Certification
- Production Capacity & Quality Control
- The actual content of the part relates only to its associated enterprise & Working Environment
- Energy Saving and Emission Reduction
- Photo

General Comments:
Chongqing Iyascare Medical Co., Ltd. was founded in 2010, is a trading company with 10 employees, the company is located at Room 406, Xiyang Free Trade Building, Shapingba District, Chongqing City, China. The company mainly supply Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products. The company has successful foreign trading experience in North America, South America etc.

Signature: [Signature] Date: 16 Sep 2021

Handcycle Wheelchair Attachment



YSH1001



FEATURES:

- 1.LED display Four different levels:
Level 1:10KM/h
Level 2:16KM/h
Level 3:20KM/h
Level 1:6KM/h (Back)
- 2.Snow leopard connection method, suitable for all normal wheelchairs (18MM-33MM)
- 3.Front light (completely waterproof)

PARAMETERS:

Driving Range:15KM
Max Speed:8KM/h
Battery: 20Ah lithium
Cushion: Artificial leather
Climb speed: 15 stair steps/Min
Climbing ability: 500 stairs steps
Stair width required: >700mm
Stair depth required: >1100mm
Stair slope required:<45°

YSH1002



PARAMETERS:

Battery indicator
USB charger
LED headlight
Product size: 450*200*180mm
Motor: 36v 180w/ 36v 250W brushless hub motor
Battery: 5Ah/ 10Ah/ 12Ah lithium
Range:15-30KM
Speed:
Level 1: 3km/h
Level 2: 6km/h
Level 1: 15km/h
N/G: 8kg/ 10kg
Packing: 570*420*320mm

YSH1003



FEATURES:

Motor: 24V 300W brushless motor
Battery: 8 Ah lithium
Charge time: 2-4 H
Brake: Drum brake
Product size: 96*48*88 cm
Range:20-30KM
Wheel size: 8 inches wheel
Frame material: aluminum
N/G:16kg/19kg
Packing: 112*27*28 cm

PARAMETERS:

- 1.LCD display different levels:
Level 1:6KM/h
Level 2:15KM/h
Level 3:20KM/h
- 2.Connection method, suitable for most manual wheelchairs (18MM-33MM)
3. LED headlight (completely waterproof)

YSH1004



FEATURES:

- 1.LED display Four different levels:
Level 1:10KM/h
Level 2:16KM/h
Level 3:20KM/h
Level 1:6KM/h (Back)
- 2.Snow leopard connection method, suitable for all normal wheelchairs (18MM-33MM)
- 3.Front light (completely waterproof)

PARAMETERS:

Wheel size: 16" 1/2 *2.4 wheels
Motor: 500W
Battery: 48V 15Ah lithium
Parking Brake: Double disc brake
Range:25-40KM
Max speed: 20KM/h
N/G weight: 14.5KG/16KG
50*25*17cm 4.85kg (charger and battery)
75*31*37cm 15.1kg
Wooden box 82*38*46cm 34kg

Handcycle Wheelchair Attachment

YSH1006



FEATURES:

- 1.LCD display Four different levels:
Level 1: 7KM/h
Level 2: 19KM/h
Level 3: 28KM/h
Level 4: 6KM/h (Back)
- 2.Easy connection method, suitable for all normal wheelchairs (18MM-33MM)
- 3.Functions: constant cruise, power display, voltage display, three speed regulation, led headlights and high decibel electronic horn.

PARAMETERS:

Product size: 1050*560*560mm
Motor: 36V 400w/ 48v 600w hub brushless motor
Battery: 36V 15Ah/ 48v 12Ah lithium
Brake: Disc brake
Range:40-55KM
Max speed: 28KM/h
Wheel size: 12" 1/2*2.4 wheels
N/G:14kg/ 17kg
Packing: 1120mm*420mm*320mm
Color: White, Black, Red

YSH1007



FEATURES:

- 1.LCD display different levels:
Level 1: 6KM/h
Level 2: 18KM/h
Level 3: 25KM/h
- 2.Connection method, suitable for most manual wheelchairs (18MM-33MM)
3. LED headlight (completely waterproof)

PARAMETERS:

Motor: 36V 350W brushless hub motor
Battery: 13 Ah lithium
Charge time: 2-4 H
Brake: Disc brake and hand brake
Product size: 96*48*88 cm
Range: 30-40KM
Wheel size: 12 inches wheel
Frame material: aluminum
N/G:21kg/24kg
Packing: 115*47*32cm

YSH1008



FEATURES:

- 1.LCD display different levels:
Level 1:15KM/h
Level 2:25KM/h
Level 3:30KM/h
- 2.Connection method, suitable for most manual wheelchairs (18MM-33MM)
3. LED headlight (completely waterproof)

PARAMETERS:

Motor: 48V 800W hub motor
Battery: 12 Ah lithium
Charge time: 4-6 H
Brake: disc brake and hand brake
Product size: 96*48*88 cm
Range:40-50KM
Wheel size: 16 inches wheel
Frame material: aluminum
N/G:26kg/29kg
Packing: 120*52*32cm

YSH1009



FEATURES:

- 1.LCD display different levels:
Level 1: 10KM/h
Level 2: 17KM/h
Level 3: 25KM/h
- 2.Connection method, suitable for most manual wheelchairs (18MM-33MM)
3. LED headlight (completely waterproof)

PARAMETERS:

Motor: 36V 400W brushless motor
Battery: 16 Ah lithium
Charge time: 4-6 H
Brake: Double disc brake
Product size: 96*48*88 cm
Range:30-40KM
Max speed: 25KM/h
Wheel size: 16 inches wheel pneumatic tire
Frame material: aluminum
N/G:18kg/20kg
Packing: 108*30*62 cm