

IYASOCARE®



FOCUS
RESEARCH
DEVELOPMENT
PRODUCTION

📍 Chongqing Iyasocare Medical Co.,Ltd.
📍 Free Trade Tower, Xiyong street. Chongqing, China
✉️ iyaso@iyasocare.com
☎️ 0086-136 1835 0425


🌐 www.iyasocare.com

CHONGQING IYASOCARE MEDICAL CO.,LTD.

PRODUCT CATALOG BROCHURE

COMPANY INTRODUCTION

Chongqing IYASOCARE MEDICAL CO.,LTD. was established in 2006 and is headquartered in the beautiful mountain city of Chongqing in China. It is an integrated industry and trade company specializing in medical device R&D, production, sales and consulting services, with business all over the world. The company's main products include wheelchairs, hospital beds, electric scooters, walkers, commode chair, etc.; the company is committed to providing customers with reasonable, complete, reliable, and high-quality medical procurement programs; since the company was founded, it has served the United States, the United Kingdom, and Australia, France, Saudi Arabia, Egypt, Pakistan, Japan, South Korea and other global customers; and has obtained medical device registration certificate, medical device production license, medical device business license, ISO13485, CE, TUV, MDR, CFS and other qualifications Certificate; is a supplier of continuous strategic cooperation with China Medical Device Association.



AUTHENTICATION CERTIFICATE

Medical Management System Certification
Award
CHONGQING IYASOCARE MEDICAL CO.,LTD.
Room 305-306, Building 3, Chongqing Railway Post Office Logistics and Storage Project, Xinyang Group 15th Floor, Yuzhou Street, Shapingba District, Chongqing

The Medical management system of the above organizations has been reviewed and complied with **ISO 13485:2016** Requirements of applicable clauses of Medical management system standards

This system is valid for the Production and sales of Medical equipment, Wheelchair, Manual wheelchair, Electric wheelchair, Scooter, Electric handcycle, Manual bed, Electric Bed, Medical bed products

Under the condition that the certificate holder's management system continues to meet the management system standards, the certification deadline is three years. And query on www.ccc-iso.com. The validity of this certificate must be confirmed by the center through regular supervision and review.

Issue date: December 26, 2023 Valid period: November 25, 2028

Logos: IAF, UKAS, France, Spain European Certified organization Limited, ANAB, IVE, AP.

DECLARATION OF CONFORMITY
ACCORDING TO (EU) 2017/745 MEDICAL DEVICE REGULATION

Manufacturer
Name: Chongqing Iyosocare Medical Co., Ltd
Address: 406, 4/F, Software Park Phase 2, 8 Xiyuan North Street (Free Trade Tower), The Netherlands
SRN: NL-AR-00000247

Product Information
Name: ELECTRIC WHEELCHAIR
Model: see annex
EMDN: Y122127
Basic UDI-DI: /
Classification: Class I, According to Rule 13, Annex VII, Regulation (EU) 2017/745

Declaration
We herewith declare that the above-mentioned products meet the requirements of Medical Device Regulation (EU) 2017/745 and the applicable standards above.

Signature: [Signature] Date: 2023.12.26
Position: GM Place: Chongqing/China

DECLARATION OF CONFORMITY
ACCORDING TO (EU) 2017/745 MEDICAL DEVICE REGULATION

Manufacturer
Name: Chongqing Iyosocare Medical Co., Ltd
Address: 406, 4/F, Software Park Phase 2, 8 Xiyuan North Street (Free Trade Tower), The Netherlands
SRN: NL-AR-00000247

Product Information
Name: Electric scooter
Model: see annex
EMDN: Y122124
Basic UDI: /
Classification: Class I, According to Rule 13, Annex VII, Regulation (EU) 2017/745

Declaration
We herewith declare that the above-mentioned products meet the requirements of Medical Device Regulation (EU) 2017/745 and the applicable standards above.

Signature: [Signature] Date: 2023.12.28
Position: GM Place: Chongqing/China

Certificate of Conformity

Certificate No.: ECT2021P095388R-1
Applicant: CHONGQING IYASOCARE MEDICAL CO.,LTD.
Room 201-205,115 Longan Avenue, Tuzhu Town, Shapingba District, Chongqing Province, China
Product Name: Wheelchair, Manual wheelchair, Electric wheelchair, Scooter, Electric handcycle, Commode chair, Walking aids, Transfer wheelchair

Trade Name: N/A
Model: (see the following annex)
Classification: Class I (according to the Manufacturer's declaration)
Review goal: Verification of the presence of Technical Documentation compatible with the Medical Devices Regulation (EU) 2017/745

This document has been issued on a voluntary basis and upon request of the manufacturer. It is our opinion that the Technical Documentation shared with us by the manufacturer is compatible with the European Standard for Medical Devices Regulation.

The manufacturer is responsible for the CE Marking process and not exempted to carry out all necessary compliance activities. This document has been issued on the basis of the regulation on ECT Voluntary Mark for the certification of products.

Shenzhen ECT Testing Technology Co., Ltd.
3th Floor, EverBuilding Business Center, Building A, No. 26 Gonghe Industrial Road, Shenzhen, Shenzhen City, China.
Tel: +86-755-87587803 Fax: +86-755-87587803 <http://www.ectlab.com>

CE Marking

TÜVRheinland Made-in-China

Report Number: MIC-AS129874 Audit Date: 27 Aug. 2022

This report is issued by Focus Technology Co., Ltd (Made-in-China.com) and the supervising inspectorate (TUV Rheinland) to confirm that:

Company Name: Chongqing Iyosocare Medical Co., Ltd.
Address: Room 201-205, 115 Longan Avenue, Tuzhu Town, Shapingba District, Chongqing City, China
Product: Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products

has been on site audited for the following scope of activity:

- General Information
- Foreign Trade Capacity
- Product Research & Development Capacity
- Management System and Product Certification
- Production Capacity & Quality Control
- The actual content of the part relates only to its associated enterprise & Working Environment
- Energy Saving and Emission Reduction & Safety

General Comments:
Chongqing Iyosocare Medical Co., Ltd. was founded in 2010, is a trading company with 10 employees, the company is located at Room 201-205, Xinyang Free Trade Building, Shapingba District, Chongqing City, China. The company mainly supply Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products. The company has successful foreign trading experience in North America, South America etc.

Signature: [Signature] Date: 16 Feb 2024

CEIG Ministerie van Volksgezondheid, Welzijn en Sport

Report Number: MIC-AS129874 Audit Date: 27 Aug. 2022

This report is issued by Focus Technology Co., Ltd (Made-in-China.com) and the supervising inspectorate (TUV Rheinland) to confirm that:

Company Name: Chongqing Iyosocare Medical Co., Ltd.
Address: Room 201-205, 115 Longan Avenue, Tuzhu Town, Shapingba District, Chongqing City, China
Product: Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products

has been on site audited for the following scope of activity:

- General Information
- Foreign Trade Capacity
- Product Research & Development Capacity
- Management System and Product Certification
- Production Capacity & Quality Control
- The actual content of the part relates only to its associated enterprise & Working Environment
- Energy Saving and Emission Reduction & Safety

General Comments:
Chongqing Iyosocare Medical Co., Ltd. was founded in 2010, is a trading company with 10 employees, the company is located at Room 201-205, Xinyang Free Trade Building, Shapingba District, Chongqing City, China. The company mainly supply Wheelchair, Hospital Bed, Scooter, Walking Aid, Medical Protection Products. The company has successful foreign trading experience in North America, South America etc.

Signature: [Signature] Date: 16 Feb 2024

Electric Wheelchair

>>> YSE123

Classic model



FEATURES:

1. Lightweight folding design, easy to carry and travel
2. Switch to manual model(3s) or manual switch
3. Detachable footrest
4. Drop back handle
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lead Acid
 Cushion: Black/Red/Blue

>>> YSE111



>>> YSE113



>>> YSE112



>>> YSE122

FDA

SFDA



>>> YSE126



>>> YSE131



MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
44.5KG/49KG	100KG	118CM	63CM	93CM	45CM	45CM	10 inches	16 inches	39CM	50CM	72*38*74CM	91*45*84CM

Electric Wheelchair

YSE109

Manual reclining



FEATURES:

1. Switch to manual model
2. Easy to fold
3. Reclining backrest 90°-150°
4. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15 km
 Max Speed: 6km/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lead Acid
 Parking Brake: Electronic parking brake
 Cushion: Red/Black/Blue

YSE108

Electric reclining



YSE118

FDA

SFDA

Electric reclining



MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
50KG/56KG	130KG	110CM	64CM	125CM	46CM	45CM	10 inches	16 inches	75CM	54CM	79*39*74CM	91*45*84CM

YSE130M

Aluminum



FEATURES:

1. It can be carried on a plane
2. Controller charging
3. Lightweight folding design, easy to carry and travel, easy manual folding
4. Switch to manual model
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15-25KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lithium battery
 Parking Brake: electromagnetic brake
 Cushion: Black

YSE130X

Auto folding



YSE130XB

Electric reclining



MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
29KG/ 33KG	120KG	98CM	63CM	110CM	48CM	45CM	8 inches	16 inches	43CM	50CM	78*48*63CM	67*50*80CM

Electric Wheelchair

YSE110

Heated seats

Ventilated seats



FEATURES:

1. Seat and backrest can be massaged, heated and ventilated, three-level adjustment
2. Headrest with Bluetooth audio connection
3. Aluminum alloy lightweight folding frame
4. Foldable armrest and footrest.

PARAMETERS:

Charge time: 4-6H
 Driving Range:15KM
 Max Speed:4.5KM/h
 Climbing ability:12°
 Motor:180*2 brushless motor
 Battery:24V 11Ah Lithium
 Parking Brake: electromagnetic braking
 Cushion:Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
25/28KG	100KG	100CM	62CM	96CM	45CM	45CM	8 inches	8 inches	44CM	50CM	67*41*62CM	65.5*44*74CM

YSE121

FDA

SFDA



FEATURES:

Charge time: 4-6H
 Driving Range: 10-20KM
 Max Speed: 6KM/h
 Climbing ability:12°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lithium battery
 Parking Brake:electromagnetic brake
 Cushion: Black

PARAMETERS:

- 1.It can be carried on a plane
- 2.Controller charging
- 3.Lightweight folding design, easy to carry and travel ,easy manual folding
- 4.Switch to manual model
- 5.Flip-up armrest

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
29KG/33KG	120KG	95CM	63CM	98 CM	45CM	45CM	8 inches	12 inches	43CM	55CM	78*63*48CM	87*66*44 CM

YSE120

Aluminum

Carbon fiber



FEATURES:

Charge time: 4-6H
 Driving Range: 10-20KM
 Max Speed: 6KM/h
 Climbing ability:9°
 Motor: 24V 180W*2Pcs brushless motor
 Battery: 24V 12Ah lithium battery
 Parking Brake:electromagnetic brake
 Cushion: Black

PARAMETERS:

- 1.It can be carried on a plane
- 2.Controller charging
- 3.Lightweight folding design, easy to carry and travel ,easy manual folding
- 4.Switch to manual model
- 5.Flip-up armrest

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
20K/22KG	100KG	95CM	63CM	98 CM	44CM	45CM	7 inches	11 inches	/	45CM	78*63*48CM	63*40*83.5CM

YSE124

Aluminum



FEATURES:

- 1.Lightweight folding design,easy to carry and travel
- 2.Easy to portable
- 3.Detachable footrest
- 4.Drop back handle
- 5.manual or electric model switch

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15-25KM
 Max Speed: 6KM/h
 Climbing ability:15°
 Motor: 250W*2Pcs/ 150W*2 brushless motor
 Battery: 24V 12Ah lithium battery
 Parking Brake: electric parking brake
 Cushion:Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
20KG/24KG	100KG	95CM	58CM	89CM	43CM	50CM	7 inches	10 inches	39CM	50CM	72*38*74CM	62*35*75CM

Electric Wheelchair

YSE128

Steel



FEATURES:

1. Lightweight folding design, easy to carry and travel
2. Switch to manual model
3. Detachable footrest
4. Drop back handle
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15 KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 12V 12Ah (6*2 lead Acid)
 Parking Brake: electromagnetic brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
26/33KG	100KG	92CM	63CM	85CM	46CM	43CM	8 inches	12 inches	50CM	50CM	70*36*71CM	71*43*76CM

YSE134



FEATURES:

1. Universal controller
2. Electric strong power offroad wheels wheelchair for disabled and elderly people.
3. Anti-slip design, the wheelchair will not slide backwards on the slope
4. Modular design, the whole disassembly, can put the auto trunk
5. Safe low-speed start design to avoid potential safety hazards caused by excessive acceleration
6. The seat can go up 45 degrees, and the back can lie down 160 degrees with electric adjustment.
7. The foot rest can go up 70 degrees adjustable.
8. Seat rotation design, 360° seat free adjustment.

PARAMETERS:

Charge time: 8-10H
 Driving Range: 20-40 KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Battery: 24V 15Ah/30Ah lithium
 Brake type: Electromagnetic brake
 Braking distance ≤800mm
 Obstacle crossing height ≤80mm
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
60/70KG	100KG	115CM	58CM	100CM	45CM	42CM	8 inches	12 inches	/	44CM	/	/

YSE130L

Aluminum



FEATURES:

1. It can be carried on a plane
2. Controller with USB mobile phone charging
3. Lightweight folding design, easy to carry and travel, easy manual folding
4. Switch to manual model
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range: 15-25KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lithium battery
 Parking Brake: electromagnetic brake
 Cushion: Black(8CM thickened)

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
24.5KG/27.5KG	100KG	94CM	62CM	89CM	44CM	45CM	8 inches	12 inches	43CM	50CM	75*62*34CM	70*43*79CM

YSE149



FEATURES:

1. 10" McNamee wheels 13" inflatable rubber rear wheel
2. Flying wing shock absorption system
3. Automotive grade seat cushion, ergonomic design, high-end leather
4. Intelligent programmed brush controller
5. Fully automatic electromagnetic brake

PARAMETERS:

Charge time: 4-6H
 Driving Range: 15-25km
 Max Speed: 6km/h
 Climbing ability: 10°
 Battery: 24V 12/ 20A Lithium battery
 Motor: 250W x2
 Parking Brake: Electromagnetic brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
50/59KG	130KG	94.5CM	64CM	88CM	47CM	41CM	10 inches	13 inches	/	55CM	94.5*64*66CM	97*67*45 CM

Electric Wheelchair

YSE150



FEATURES:

Seat width:40 cm
 Seat depth:44cm
 Seat height:54 cm
 Product size: 102*55*99 cm
 Rear wheel size: 9 inches
 Front wheel size: 8 inches
 N/G: 30/ 35KG
 Frame material: steel
 Load:130 KG
 Packing size:

PARAMETERS:

1. 8-inch McNamee front wheel
2. aluminum alloy frame
3. Seat height adjustable, telescopic footrest, foldable armrests
4. Intelligent controller with lighting function
5. Fully automatic electromagnetic brake
6. 3D butterfly wing cushion backrest

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
30/35KG	130KG	102CM	55CM	99CM	40CM	44CM	8 inches	9 inches	/	54CM	/	/

YSE154

Steel



FEATURES:

- 1, Folding front and back
- 2, Armrests can be turned over
- 3, Pedal can be turned over
- 4, Battery quick disassembly
- 5, Controller front and back adjustment, left and right interchangeable

PARAMETERS:

Charge time: 4-6H
 Driving Range: 15KM
 Max Speed: 7KM/h
 Climbing ability:12°
 Motor: 250W*2Pcs
 Battery: 24V 12Ah lead Acid
 Parking Brake>manual brake and electric parking brake
 Cushion:Black/Red/Blue

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
33/36kg	100KG	102CM	67CM	89M	43CM	44CM	8 inches	12 inches	78CM	51CM	70*42*73CM	88*42*86CM

IYASOCARE®

Off Road Wheelchair

YSE101A

All terrain



FEATURES:

1. Upgrade the thickened and widened seat, made of top grade leather.
2. Upgrade the new thickened leather armrests.
3. The body volume is narrower and shorter.
4. The rear wheel hub has increased the night driving lights and turn signals, more secure.
5. Foot rest can be adjusted forward and backward.

PARAMETERS:

Charge time: 6-8H
 Driving Range: 25-60km
 Max Speed: 10 km/h
 Climbing ability: 15°
 Motor: 24V 250W*2
 Battery: 24V 40Ah lead acid

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
80/85KG	200KG	110CM	62.5CM	121CM	45.5CM	45CM	10 inches	14 inches	/	56CM	/	115*65*60CM

YSE103

All terrain

Double seat



FEATURES:

1. Increase the rear seat, can be traveled by two people.
2. Dual handles are easy to control, the tires are enlarged and widened, equipped with spring shock absorption.
3. Large motor power, fast speed, large battery capacity, long range.
4. Suitable for paraplegic, hemiplegic patients, lower limb disabled and mobility impaired elderly.
5. Suitable for long-distance traveling, demanding speed and dealing with poor road surface.

PARAMETERS:

Charging time: 6-8h
 Driving range : 25-75km
 Max speed: 20 km/h
 Climbing capacity : 12°
 Motor: 800W differential motor
 Battery: 48v 20ah / 32ah lead acid / 60v 30ah lithium battery
 Parking brake: drum brake
 Cushion: Red/Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
126/135KG	200KG	166CM	90CM	125CM	46CM	45CM	13 inches	14 inches	/	57CM	/	145*80*70CM

YSE102A

All terrain



FEATURES:

1. Upgraded 800w motor, max speed 25km/hour.
2. Upgraded disc brakes on both front and rear, and added oil brakes on the rear wheel.
3. Upgraded wider and thicker leather seats than the first generation.
4. Upgraded the front and rear wheels.
5. Upgraded the center integrated pedal.

PARAMETERS:

Charging time: 6-8h
 Driving range: 25-75km
 Max speed: 25 km/h
 Climbing capacity: 15°
 Motor: 800W differential motor
 Battery: 48v 20ah / 32ah lead acid / 60v 30ah lithium battery
 Parking brake: disc brake, oil brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
105/110KG	200KG	126CM	65.5CM	125CM	46CM	45CM	14 inches	16 inches	/	57CM	/	115*65*64CM

YSE104

All terrain

Electric reclining



FEATURES:

1. Universal controller
2. Electric strong power offroad wheels wheelchair for disabled and elderly people.
3. Can climb 10cm high stairs.
4. Anti-slip design, the wheelchair will not slide backwards on the slope
5. Standard USB charging port, safety belt.

PARAMETERS:

Product size: 115*68*110 cm
 Rear wheel size: 14 inches
 Front wheel size: 10 inches
 Frame material: steel
 Load: 200KG
 Packing size: 115*65*60cm
 Wooden box: 120*72*79cm 140KG

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
70/80KG	200KG	115CM	68CM	110CM	45 CM	45CM	10 inches	14 inches	78CM	52CM	/	115*65*60CM

Off Road Wheelchair

YSE105

All terrain



FEATURES:

1. Electric strong power offroad wheels wheelchair for disabled and elderly people.
2. Standard USB charging port, safety belt, optional selfie stick bracket and personalized storage basket.
3. Speed and battery life will vary with the battery & motor power changes.
4. Four driving modes are available: off-road, climb, normal, and indoor.
5. 180° reclining adjustable car seat.

PARAMETERS:

Charging time: 4-6h
 Driving range :40-60km
 Top speed :28 km/h
 Battery: 48v 20ah lead acid
 Parking brake: disc brake

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	200KG	100CM	63CM	58CM	49CM	46CM	14 inches	14inches	/	58CM	/	135*70*90CM

YSE107

All terrain



FEATURES:

1. Suitable for disabled and elderly people's electric powerful cross-country wheelchair.
2. It can climb stairs of 10 centimeters high.
3. Anti-slip slope design, the wheelchair will not slide backward or roll over on the slope.
4. Comes standard with USB charging port, leather seat, seat belt, turn signal, headlight and horn.
5. Speed and battery life vary with battery and motor power.
6. Four driving modes are available: off-road, climb, normal, and indoor.

PARAMETERS:

Charging time: 4-6h
 Driving range :40-60km
 Top speed :28 km/h
 Climbing capacity :15°
 Battery: 48v 20ah lead acid
 Parking brake: disc brake
 Cushion:black

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	200KG	125CM	64CM	108 CM	49CM	46CM	14 inches	14 inches	/	58CM	/	135*70*90CM

IYASOCARE®

Small Light Electric Wheelchair

YSE119

Aluminum



FEATURES:

1. High-strength aluminum alloy frame, sturdy and durable
2. Brushless electromagnetic brake motor, safe and no slippery slope
3. Low noise Li-ion battery, lightweight and convenient
4. Long life brushless motor and controller, 360 degree flexible control.
5. Flip-up armrest

PARAMETERS:

Charge time: 4-6H
 Driving Range:15 KM
 Max Speed:6KM/h
 Climbing ability:8°
 Motor: 180w*2 brushless motor
 Battery:24V 6Ah Lithium
 Parking Brake>manual brake and electronic parking brake
 Cushion:Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
20/22KG	100KG	112CM	60CM	94CM	43CM	45CM	7 inches	12 inches	44CM	48CM	40*60*84CM	43*62*87CM

YSE136

Ergonomic seat



FEATURES:

Charge time: 4-6H
 Driving Range:5-10KM
 Max Speed:6KM/h
 Climbing ability:9°
 Motor:250W*2Pcs
 Battery:24V 6Ah/ 10Ah/ 20Ah lithium battery
 Parking Brake:electromagnetic brake
 Cushion: Black

PARAMETERS:

1. Aluminum alloy frame, sand black finish.
2. Ergonomic seat plate with removable seat cushion.
3. Front wheel shock absorber.
4. Rear flip armrest with PU armrest pad.
5. Adjustable leg height.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
20KG/22KG	100KG	105CM	59CM	90CM	41CM	40CM	7 inches	10 inches	/	48CM	/	61*36*80CM

YSE132

Shock absorbers



FEATURES:

Charge time: 4-6H
 Driving Range:5-10KM
 Max Speed:6KM/h
 Climbing ability:9°
 Motor:250W*2Pcs
 Battery:24V 6Ah/ 10Ah/ 20Ah lithium battery
 Parking Brake:electromagnetic brake
 Cushion: Black

PARAMETERS:

1. Aluminum alloy frame.
2. Breathable backrest, high quality leather padding.
3. Front wheel shock absorber
4. Rear flip armrest is equipped with PU armrest pad.
5. Front and rear wheels are equipped with PU non-inflatable honeycomb wheels.
6. Adjustable leg height

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
22KG/24KG	100KG	117CM	53CM	93M	41CM	40CM	8 inches	10 inches	/	48CM	/	55*34*82 CM

YSE142

Magnesium alloy



FEATURES:

1. It can be carried on a plane
2. Controller charging
3. Lightweight folding design, easy to carry and travel ,easy manual folding
4. Switch to manual model
5. Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range:15 KM
 Max Speed: 6KM/h
 Climbing ability: 9°
 Motor: 150W*2Pcs brushless motor
 Battery:24V 6Ah/12Ah lithium
 Cushion:Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
18.5/22.5KG	100KG	92CM	60CM	93M	43CM	44CM	7 inches	10 inches	48CM	51CM	60*28*70CM	62*35*75CM

Small Light Electric Wheelchair

YSE143

Ergonomic seat



FEATURES:

Charge time: 4-6H
 Driving Range: 5-10KM
 Max Speed: 6KM/h
 Climbing ability: 9°
 Motor: 250W*2Pcs brushless motor
 Battery: 24V 6Ah/ 10Ah/ 20Ah lithium battery
 Parking Brake: electromagnetic brake
 Cushion: Black

PARAMETERS:

1. Spraying sand black + spraying orange color
2. Comfortable ergonomic seat plate with premium set of removable leather cushions
3. 8 inch front wheel with PU+shock absorber
4. Rear flip armrest with PU armrest pads
5. Leg height adjustable

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
18.5 KG/22 KG	100KG	115CM	59CM	87M	41CM	40CM	8 inches	10 inches	/	48CM	/	61*36*80CM

YSE147

Ergonomic seat



FEATURES:

Charge time: 4-6H
 Driving Range: 5-20KM
 Max Speed: 6KM/h
 Climbing ability: 9°
 Motor: 180W*2Pcs
 Battery: 24V 6Ah/ 10Ah/ 20Ah lithium battery
 Parking Brake: electromagnetic brake
 Cushion: Black

PARAMETERS:

1. Spraying sand black
2. Ergonomic second generation backrest, seat plate with advanced set of removable leather cushion.
3. 8 inch front wheel with PU+shock absorber.
4. 12 inch rear explosive-proof honeycomb run-flat tire
5. Rear flip armrest is equipped with PU armrest pad.
6. Configuration headrest

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
23KG/25KG	100KG	112CM	59CM	90M	40CM	40CM	8 inches	12 inches	50CM	48CM	97*59**60CM	61*36*80CM

YSE143A

Ergonomic seat



FEATURES:

Charge time: 4-6H
 Driving Range: 5-10KM
 Max Speed: 6KM/h
 Climbing ability: 9°
 Motor: 250W*2Pcs
 Battery: 24V 6Ah/ 10Ah/ 20Ah lithium battery
 Parking Brake: electromagnetic brake
 Cushion: Black

PARAMETERS:

1. Spraying sand black
2. Ergonomic second generation backrest, seat plate with advanced set of removable leather cushion.
3. 8 inch front wheel with PU+shock absorber.
4. Rear flip armrest is equipped with PU armrest pad.
5. Adjustable leg height
6. Configuration headrest

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
18.5 KG/22 KG	100KG	115CM	59CM	87M	41CM	40CM	8 inches	10 inches	/	48CM	/	61*36*80CM

IYASOCARE®

Carbon Fiber Electric Wheelchair

YSE139

Carbon fiber



FEATURES:

- 1.It can be carried on a plane
- 2.Lightweight folding design, easy to carry and travel
- 3.Carbon fiber frame material
- 4.Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range:10-15 KM
 Max Speed:4.5KM/h
 Climbing ability: $\geq 3^\circ$
 Motor: 150W*2Pcs
 Battery:24V 6Ah/10Ah
 Braking distance ≤ 1500 mm
 Obstacle crossing height ≤ 25 mm
 Cushion:Black

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	115CM	96CM	86CM	40CM	40CM	6 inches	12inches	24CM	43CM	74*36*86CM	89*42*80CM

YSE140

Carbon fiber



FEATURES:

- 1.It can be carried on a plane
- 2.Lightweight folding design, easy to carry and travel
- 3.Carbon fiber frame material
- 4.Flip-up armrest

PARAMETERS:

Charge time: 6-8H
 Driving Range:15 KM
 Max Speed: 6KM/h
 Climbing ability: 9°
 Motor: 220W*2Pcs brushless motor
 Battery:24V 6Ah/12Ah lithium
 Cushion:Black

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	93CM	59CM	92CM	43CM	41CM	7 inches	10 inches	49CM	50CM	59*29*73CM	62*35*75CM

IYASOCARE®

Electric Stair Climbing Wheelchair

YSE201



FEATURES:

1. Universal controller
2. High quality stair climber stair climbing chair for disabled and elderly people
3. Climbing lock: slow down and reverse lock
4. Flat ground and stair climbing dual use

PARAMETERS:

Driving Range: 15KM
 Max Speed: 8KM/h
 Battery: 20Ah lithium
 Cushion: Artificial leather
 Climb speed: 15 stair steps/Min
 Climbing ability: 500 stairs steps
 Stair width required: >700mm
 Stair depth required: >1100mm
 Stair slope required: <45°

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	100CM	65CM	120CM	44CM	45CM	8 inches	12 inches	53CM	50CM	98*65*50CM	112*75*79CM

YSE203



FEATURES:

1. Universal controller
2. Multifunctional height adjustable folding wheelchair for stair climbing
3. Textile breathable sponge cushion, pillow adjustable
4. Electric driving and electric climbing dual use

PARAMETERS:

Driving Range: 15KM
 Max Speed: 8KM/h
 Battery: 20Ah lithium
 Driving motor: 24V 250W*2
 Climbing motor: 24V 250W
 Cushion: Artificial leather
 Climb speed: 15 stair steps/Min
 Climbing ability: 500 stairs steps
 Stair width required: >700mm
 Stair depth required: >1100mm
 Stair slope required: <45°

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	102CM	70CM	110CM	44CM	44CM	7.5 inches	12 inches	78CM	44CM	70*58*93CM	112*75*79CM

YSE202



FEATURES:

Battery: 24V 13Ah lithium
 Cushion: Artificial leather
 Climb speed: 30 stair steps/Min
 Climbing ability: 500 stairs steps
 Stair width required: >650mm
 Staircase corner width and depth required: >1200mm
 Stair slope required: <50°

PARAMETERS:

1. Aluminum alloy material, light weight, can be folded
2. Wide range of applications, can load most of the manual wheelchairs on the market
3. Convenient operation, stable climbing, no high threshold for use
4. small size, suitable for most residential stairs

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	160CM	77CM	70CM	/	/	2.5 inches	5 inches	/	/	105*50*23 CM	130*56*18CM

YSE204



FEATURES:

1. Lightweight multifunctional height adjustable folding wheelchair for stair climbing
2. 3 gears climbing speed
3. New large wheel, easy to deal with various road conditions
4. Four-season dual-use seat cushion

PARAMETERS:

Battery: 13 Ah lithium
 Cushion: Artificial leather
 Climb speed: 15 stair steps/Min
 Climbing ability: 500 stairs steps
 Stair width required: >650mm
 Stair depth required: >1100mm
 Stair slope required: <45°

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	108CM	50CM	151CM	45CM	45CM	4 inches	12 inches	50CM	50CM	50*29*103CM	104*54*30CM

Electric Stair Climbing Wheelchair

YSE205



FEATURES:

1. Universal controller
2. High quality stair climber stair climbing chair for disabled and elderly people
3. Climbing lock: slow down and reverse lock
4. Flat ground and stair climbing dual use
5. Reclining chair

PARAMETERS:

Driving Range: 15KM
 Max Speed: 8KM/h
 Motor: 24V DC 370W
 Battery: 20Ah (lithium battery)
 Parking brake: manual brake and electromagnetic brake
 Cushion: Artificial leather
 Minimum depth of stairs: 1.3m
 Climbed Number: 500
 Stair width: $\geq 0.7m$
 Stair slope: $< 45^\circ$
 Climb speed: 15 stair steps/Min

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
93/100KG	150KG	108CM	68CM	130CM	48CM	43CM	8 inches	12 inches	/	44CM	68*77*88CM	/

YSE208



FEATURES:

1. Universal controller
2. Multifunctional height adjustable folding wheelchair for stair climbing
3. Electric driving and electric climbing dual use

PARAMETERS:

Driving Range: 15KM
 Battery: 13Ah lithium
 Cushion: Artificial leather
 Climb speed: 15 stair steps/Min
 Climbing ability: 500 stairs steps
 Stair width required: $> 700mm$
 Stair depth required: $> 1100mm$
 Stair slope required: $< 45^\circ$

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
38KG/43KG	150KG	98CM	65CM	125CM	43CM	45CM	8 inches	12 inches	/	44CM	110*65*35CM	123*60*71CM

YSE207



FEATURES:

1. Universal controller
2. High quality stair climber stair climbing chair for disabled and elderly people
3. Climbing lock: slow down and reverse lock
4. Flat ground and stair climbing dual use
5. Backup camera
6. Climb the stairs independently
7. Large comfortable seat

PARAMETERS:

Driving Range: 15KM
 Battery: 20Ah lithium
 Cushion: Artificial leather
 Climb speed: 15 stair steps/Min
 Climbing ability: 500 stairs steps
 Stair width required: $> 850mm$
 Stair depth required: $> 1500mm$
 Stair slope required: $< 45^\circ$

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
130KG	150KG	128CM	76CM	140CM	50CM	50CM	8 inches	12 inches	/	63CM	114*76*120CM	78*126*125CM

YSE209



FEATURES:

- Battery: 10Ah lithium
 Cushion: Artificial leather
 Climb speed: 15 stair steps/Min
 Climbing ability: 500 stairs steps
 Stair width required: $> 700mm$
 Staircase corner width and depth required: $> 1100mm$
 Stair slope required: $< 45^\circ$

PARAMETERS:

1. Electrical stair climbing
2. Textile breathable sponge cushion
3. Manual drive and electrical stair climbing dual use

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
27KG/31KG	100KG	81CM	57CM	145CM	41CM	35CM	4 inches	5 inches	41CM	/	102*57*34CM	104*58*35CM

Electric Stair Climbing Wheelchair



YSE-F1

FEATURES:














1. Simple operation, can be quickly learned to use.
2. It can be docked with manual wheelchairs of different types, with wide versatility.

PARAMETERS:

Battery: 12Ah
Climb speed: 15 stair steps/Min
Climbing ability: 900 stairs steps (1 hour)
Stair width required: >700mm
Stair depth required: >1300mm
Stair slope required: <45°



MAIN PARAMETERS

 N/G	 Load	 overall length	 overall width	 overall height	 seat width	 seat depth	 front wheel	 rear wheel	 backrest height	 seat height	 folding size	 DIM
80KG/98KG	130KG	102CM	45CM	29CM	/	/	/	/	/	44CM	98*65*50CM	112*75*79CM

IYASOCARE®

Standing Wheelchair

YSE301

Full lying



FEATURES:

1. Electromagnetic brake system, carbon steel body,
2. Hydraulic adjustment device can be completely lay down.
3. Front and rear shock absorption design, rocker anti-radio interference controller.
4. Equipped with high-quality knee pads.

PARAMETERS:

Driving Range: 45KM
 Max Speed: 8KM/h
 Battery: 24V 50A(lead acid)
 Motor: 24V 350W*2
 Parking Brake: Intelligent electronic brake + electromagnetic brake
 Cushion: Black

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	120CM	65CM	120CM	48CM	45CM	14 inches	10 inches	52CM	55CM	/	125*70*70CM

YSE304

Full lying



FEATURES:

1. Electromagnetic brake system, carbon steel body.
2. It can be electrically adjusted to standing mode and full lying mode, manual electric dual-mode system, multi-function remote control operating panel, each function button exists independently, optional toggle switch, operated by wrist.
3. Large-capacity battery, long life use.
4. Equipped with high quality knee pads and double-sided lifting armrests.

PARAMETERS:

Driving Range: 45KM
 Max Speed: 8KM/h
 Battery: 24V 55Ah
 Motor: 24V 450W
 Parking Brake: Intelligent electronic brake + electromagnetic brake
 Cushion: Black

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	120KG	120CM	65CM	110CM	45CM	45CM	8 inches	15 inches	/	55CM	/	112*72*95CM

YSE302

All terrain



FEATURES:

1. Suitable for disabled and elderly people's electric powerful cross-country wheelchair.
2. It can climb stairs of 10 centimeters high.
3. Anti-slip slope design, the wheelchair will not slide backward or roll over on the slope.
4. Comes standard with USB charging port, leather seat, seat belt, turn signal, headlight and horn.
5. Speed and battery life vary with battery and motor power.
6. Three driving modes are available: off-road, normal, and indoor.
7. The seat can be adjusted to 180 degrees of recline.
8. Motorized standing available with remote control to assist standing anytime, anywhere.

PARAMETERS:

Charging time: 8-10h
 Driving range :40-60km
 Top speed :28 km/h
 Climbing capacity :15°
 Battery: 48v 20ah lead acid
 Parking brake: electromagnetic brake
 Cushion:black

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	200KG	125CM	64CM	108CM	49CM	46CM	14 inches	14 inches	/	/	/	135*70*90CM

YSE305



FEATURES:

1. Electric walking, electric standing, electric standing and walking
2. Turn on the spot with 360 degrees
3. One-hand joystick controller
4. Universal front wheel
5. Armrests on both sides can be raised

PARAMETERS:

Driving Range:30-40 KM
 Max Speed:8KM/h
 Battery:24V 40A(lead acid)
 Motor:250W*2 24V
 Parking Brake:Electromagnetic Brake and intelligent electronic brake
 Cushion:Black with red

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	100CM	56CM	85M	45CM	43CM	8 inches	14/22 inches	/	47CM	/	115*65*60CM

Standing Wheelchair

YSE306

All terrain



FEATURES:

Charging time: 6-8h
 Driving range :25-60km
 Motor: 48V 600W motor
 Top speed :20 km/h
 Climbing capacity :15°
 Battery: 48v 20ah lead acid
 Parking brake: drum brake
 Cushion:black/ red

PARAMETERS:

1. Electric off-road wheelchair suitable for disabled and elderly people.
2. Slope anti-slip design
3. Standard equipped with standing controller
4. 3 speed adjustable, maximum speed 20km/h. 8.
5. Motorized standing by remote control

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	120KG	125CM	65CM	125M	45CM	49CM	14 inches	13 inches	/	53CM	/	125*70*68CM

YSE308

Full lying



FEATURES:

1. Electromagnetic brake system, carbon steel body.
2. 4 modes: seating/elevating/reclining/standing. Multi-function remote control operation panel, each function button is independent, optional toggle switch, operated by wrist.
3. Large capacity battery with long service life.
4. Equipped with high quality knee pads and double-sided lifting armrests.

PARAMETERS:

Driving Range: 45KM
 Max Speed: 6KM/h
 Gradeability: 12°
 Battery: 24V 55Ah
 Motor: 320w*2 worm gear motor
 Parking Brake: Intelligent electronic brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	120KG	132CM	65CM	133M	43CM	46CM	8 inches	14 inches	54CM	56-69CM	/	118*71*86CM

YSE309



FEATURES:

1. Electromagnetic brake system, carbon steel body,
2. Hydraulic adjustment device can be completely lay down.
3. Front and rear shock absorption design, rocker anti-radio interference controller.
4. Equipped with high-quality knee pads, and double-sided lifting armrests.

PARAMETERS:

Driving Range:8-15km
 Max Speed:0-8 km/h
 Battery:20A lithium battery
 Parking Brake:hand brake + electromagnetic brake
 Cushion:Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
100KG/130KG	150KG	108CM	68CM	130CM	49CM	45CM	12 inches	8 inches	/	62CM	/	/

YSE311

Semi-automatic standing



FEATURES:

- Power standing wheelchair, armrest, comfortable seat belt, power control standing, manual control of direction and progress.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
34.5/39KG	100 KG	/	65CM	104CM	46CM	/	8 inches	24 inches	/	58CM	/	110*69*67CM

Standing Wheelchair

>>> YSM4005














Manual standing



FEATURES:

1. Manual standing up
2. Comfortable cushion
3. Free rubber big rear wheel with manual ring
4. Steel frame, adjustable wheelchair height

MAIN PARAMETERS

 N/G	 Load	 overall length	 overall width	 overall height	 seat width	 seat depth	 front wheel	 rear wheel	 backrest height	 seat height	 folding size	 DIM
26KG/34KG	100KG	106CM	67CM	85CM	42CM	43CM	5 inches	24 inches	/	45CM	/	115*65*64CM

IYASOCARE®

Electric Lifting Wheelchair

YSE307



FEATURES:

1. One-hand joystick controller, front wheel universal steering
2. The armrests on both sides can be quickly removed, the backrest can be reclined, the pedal can be quickly removed, and the lifting pedal can be replaced
3. Turned on the spot in a 360-degree,
4. It has electric walking and electric lifting functions

PARAMETERS:

Driving Range: 25 km
 Max Speed: 0-8 km/h
 Motor: 24V 250W*2
 Battery: 24V 24A lead acid
 Parking Brake: Intelligent electronic brake + electromagnetic brake
 Cushion: Black and red

MAIN PARAMETERS

Icon	Parameter	Value
	N/G	N/G
	Load	150KG
	overall length	105CM
	overall width	59CM
	overall height	110-165CM
	seat width	45CM
	seat depth	45CM
	front wheel	8 inches
	rear wheel	12 inches
	backrest height	72CM
	seat height	45-95CM
	folding size	/
	DIM	/

YSE314



FEATURES:

1. One-hand joystick controller, front wheel universal steering
2. Luxurious and comfortable car seats for long periods of use
3. It has electric driving and electric lifting functions

PARAMETERS:

Charge time: 8-10H
 Driving Range: 15 KM
 Max Speed: 6KM/h
 Climbing ability: 12°
 Motor: 24V 180W*2
 Battery: 20AH lithium
 Brake type: Electromagnetic brake
 Braking distance ≤ 800mm
 Obstacle crossing height ≤ 40mm
 Cushion: Black

MAIN PARAMETERS

Icon	Parameter	Value
	N/G	N/G
	Load	125KG
	overall length	93CM
	overall width	59CM
	overall height	140CM
	seat width	47CM
	seat depth	45CM
	front wheel	8 inches
	rear wheel	12.5 inches
	backrest height	/
	seat height	/
	folding size	75*59*69CM
	DIM	77*63*72KG

IYASOCARE®

Handcycle Wheelchair Attachment

YSH1001



FEATURES:

- 1.LED display Four different levels:
Level 1:10KM/h
Level 2:16KM/h
Level 3:20KM/h
Level 1:6KM/h (Back)
- 2.Snow leopard connection method, suitable for all normal wheelchairs (18MM-33MM)
- 3.Front light (completely waterproof)

PARAMETERS:

Driving Range:15KM
Max Speed:8KM/h
Battery: 20Ah lithium
Cushion: Artificial leather
Climb speed: 15 stair steps/Min
Climbing ability: 500 stairs steps
Stair width required: >700mm
Stair depth required: >1100mm
Stair slope required:<45°

YSH1002



PARAMETERS:

Battery indicator
USB charger
LED headlight
Product size: 450*200*180mm
Motor: 36v 180w/ 36v 250W brushless hub motor
Battery: 5Ah/ 10Ah/ 12Ah lithium
Range:15-30KM
Speed:
Level 1: 3km/h
Level 2: 6km/h
Level 1: 15km/h
N/G: 8kg/ 10kg
Packing: 570*420*320mm

YSH1003



FEATURES:

Motor: 24V 300W brushless motor
Battery: 8 Ah lithium
Charge time: 2-4 H
Brake: Drum brake
Product size: 96*48*88 cm
Range:20-30KM
Wheel size: 8 inches wheel
Frame material: aluminum
N/G:16kg/19kg
Packing: 112*27*28 cm

PARAMETERS:

- 1.LCD display different levels:
Level 1:6KM/h
Level 2:15KM/h
Level 3:20KM/h
- 2.Connection method, suitable for most manual wheelchairs (18MM-33MM)
3. LED headlight (completely waterproof)

YSH1004



FEATURES:

- 1.LED display Four different levels:
Level 1:10KM/h
Level 2:16KM/h
Level 3:20KM/h
Level 1:6KM/h (Back)
- 2.Snow leopard connection method, suitable for all normal wheelchairs (18MM-33MM)
- 3.Front light (completely waterproof)

PARAMETERS:

Wheel size: 16" 1/2 *2.4 wheels
Motor: 500W
Battery: 48V 15Ah lithium
Parking Brake: Double disc brake
Range:25-40KM
Max speed: 20KM/h
N/G weight: 14.5KG/16KG
50*25*17cm 4.85kg (charger and battery)
75*31*37cm 15.1kg
Wooden box 82*38*46cm 34kg

Handcycle Wheelchair Attachment

YSH1006



FEATURES:

- 1.LCD display Four different levels:
Level 1: 7KM/h
Level 2: 19KM/h
Level 3: 28KM/h
Level 4: 6KM/h (Back)
- 2.Easy connection method, suitable for all normal wheelchairs (18MM-33MM)
- 3.Functions: constant cruise, power display, voltage display, three speed regulation, led headlights and high decibel electronic horn.

PARAMETERS:

Product size: 1050*560*560mm
Motor: 36V 400w/ 48v 600w hub brushless motor
Battery: 36V 15Ah/ 48v 12Ah lithium
Brake: Disc brake
Range:40-55KM
Max speed: 28KM/h
Wheel size: 12" 1/2*2.4 wheels
N/G:14kg/ 17kg
Packing: 1120mm*420mm*320mm
Color: White, Black, Red

YSH1007



FEATURES:

- 1.LCD display different levels:
Level 1: 6KM/h
Level 2: 18KM/h
Level 3: 25KM/h
- 2.Connection method, suitable for most manual wheelchairs (18MM-33MM)
3. LED headlight (completely waterproof)

PARAMETERS:

Motor: 36V 350W brushless hub motor
Battery: 13 Ah lithium
Charge time: 2-4 H
Brake: Disc brake and hand brake
Product size: 96*48*88 cm
Range: 30-40KM
Wheel size: 12 inches wheel
Frame material: aluminum
N/G:21kg/24kg
Packing: 115*47*32cm

YSH1008



FEATURES:

- 1.LCD display different levels:
Level 1:15KM/h
Level 2:25KM/h
Level 3:30KM/h
- 2.Connection method, suitable for most manual wheelchairs (18MM-33MM)
3. LED headlight (completely waterproof)

PARAMETERS:

Motor: 48V 800W hub motor
Battery: 12 Ah lithium
Charge time: 4-6 H
Brake: disc brake and hand brake
Product size: 96*48*88 cm
Range:40-50KM
Wheel size: 16 inches wheel
Frame material: aluminum
N/G:26kg/29kg
Packing: 120*52*32cm

YSH1009



FEATURES:

- 1.LCD display different levels:
Level 1: 10KM/h
Level 2: 17KM/h
Level 3: 25KM/h
- 2.Connection method, suitable for most manual wheelchairs (18MM-33MM)
3. LED headlight (completely waterproof)

PARAMETERS:

Motor: 36V 400W brushless motor
Battery: 16 Ah lithium
Charge time: 4-6 H
Brake: Double disc brake
Product size: 96*48*88 cm
Range:30-40KM
Max speed: 25KM/h
Wheel size: 16 inches wheel pneumatic tire
Frame material: aluminum
N/G:18kg/20kg
Packing: 108*30*62 cm

Foldable Lightweight Wheelchair

YSM2026



FEATURES:

1. Aluminum alloy frame, sandblack + carbon transfer treatment, rear dual large shock absorbers
2. Ergonomic seat back with removable waterproof cushion.
3. front wheel shock absorber
4. Linkage rear flip armrest with PU armrest pads
5. Self-locking brake
6. Front and rear wheels are equipped with honeycomb non-inflatable shock-absorbing wheels.
7. One button unlocking
8. Adjustable footrest

Aluminum

Ergonomic seat

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	88CM	55CM	92CM	41CM	40CM	6 inches	12 inches	/	50CM	45*50*73 CM	55*34*82 CM

YSM2027



FEATURES:

1. Aluminum alloy frame, sand black + carbon fiber transfer treatment.
2. Ergonomics II backrest with removable waterproof cushion on the seat back.
3. Front wheel shock absorber
4. Linkage rear flip armrest with PU armrest pads
5. Configure self-locking brake
6. front and rear wheels are equipped with honeycomb non-inflatable shock-absorbing wheels.
7. One-touch unlocking
8. Adjustable footrest
9. Hydraulic reclining adjustment

Aluminum

Ergonomic seat

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	96CM	55CM	98CM	41CM	40CM	11 inches	12 inches	/	50CM	45*50*73 CM	55*34*82 CM

YSM2028



FEATURES:

1. Aluminum alloy lightweight frame
2. Ergonomic seat backrest, wrapped in mesh black elemental fabric.
3. Aluminum alloy frame with plastic footrests that can be turned up and adjusted.
4. 6-inch wheels, shock absorbing fork. Rear 12-inch wheels
5. With rear hand brake
6. Wheelchair can be folded, armrests can be upturned and fixed, aluminum alloy footrest can be upturned.
7. Telescopic drawbar, easy to pull after folding.
8. Rear push handle is designed.

Aluminum

Ergonomic seat

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	102CM	52.5CM	88CM	42CM	45CM	6 inches	12 inches	43CM	50CM	45*50*73 CM	57*41*85 CM

YSM2029



FEATURES:

1. Aluminum alloy lightweight frame
2. Ergonomic seat backrest, wrapped in mesh black elemental fabric.
3. Aluminum alloy frame with plastic footrests that can be turned up and adjusted.
4. 6-inch honeycomb tires on front wheel, shock absorbing fork. Rear 12-inch honeycomb tires
5. With rear hand brake
6. Wheelchair can be folded, armrests can be upturned and fixed, aluminum alloy footrest can be upturned.
7. Rear push handle is designed and detachable.

Aluminum

Ergonomic seat

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N	100KG	102CM	55CM	88CM	42CM	45CM	6 inches	12 inches	/	46CM	40*55*70 CM	57*39*78 CM

Foldable Lightweight Wheelchair

YSM2030
















FEATURES:

1. Aluminum alloy frame
2. Light weight
3. Double brakes
4. Flip-up footrest. Backrest with double insurance
5. Flip-up armrest

Aluminum

MAIN PARAMETERS

 N	 Load	 overall length	 overall width	 overall height	 seat width	 seat depth	 front wheel	 rear wheel	 backrest height	 seat height	 folding size	 DIM
6.9KG	75KG	74CM	52CM	88CM	34CM	32CM	6 inches	8 inches	41CM	49CM	72*53*32CM	71*30*52CM

IYASOCARE®

Manual Wheelchair - Standard Wheelchair

YSM2001

TYPE: A



FEATURES:

1. Foldable backrest
2. Steel handrail with movable leg folding handle
3. PU handrail, high-grade camouflage 3D breathable seat cushion
4. High quality nylon pedal, protective plate color aluminum plate

Carbon Steel Frame

N/G: 12/14KG

TYPE: B



FEATURES:

1. Aluminum alloy frame light weight
2. Movable handrails, footrest
3. PU armrest, high-grade camouflage 3D breathable, seat cushion, high quality nylon pedal, color aluminum plate
4. 20" aluminum alloy PU tire, equipped with blade light.

Aluminum Frame

N/G: 11/12.8KG

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
/	100KG	95CM	60CM	85CM	44CM	37CM	6 inches	20 inches	65CM	45CM	80*24*70CM	83*26*77CM

YSM2004

TYPE: A



FEATURES:

1. Carbon steel frame, strong and durable
2. Movable handrails, footrest
3. PU armrest, high-grade camouflage 3D breathable, seat cushion, high quality nylon pedal, color aluminum plate

Carbon Steel

N/G: 13/16KG

TYPE: B



FEATURES:

1. Aluminum alloy frame light weight
2. Movable handrails, footrest
3. PU armrest, high-grade camouflage 3D breathable, seat cushion, high quality nylon pedal, color aluminum plate
4. 22" aluminum alloy wheels

Aluminum Frame

N/G: 12/14.8KG

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
/	100KG	104CM	61CM	88CM	46CM	41CM	6 inches	22 inches	/	/	81*26*86CM	82*27*88CM

Manual Wheelchair - Standard Wheelchair

YSM2005



FEATURES:

1. Carbon steel frame with gray painted finish
2. reversible armrests, PVC leather armrest pads
3. Detachable footrests, height adjustable.
4. Seat back cushion made of PVC leather, backrest can be folded.
5. 8-inch wider front casters, 24-inch aluminum rear wheels, wider mountain tires
6. Elbow joint brakes, lower than the seat surface after braking.

Fat Tire

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
20KG/22KG	140KG	100CM	72CM	94CM	56CM	41CM	8 inches	24 inches	40CM	/	100*26*94 CM	81*28*88 CM

YSM2012



FEATURES:

1. Aluminum alloy two-tone
Frame: Aluminum alloy two-color frame, with frosted sprayed surface.
2. Oxford nylon fabric with mesh cushion, seat back cushion
Seat back pad: Oxford nylon fabric with mesh cushion, seat back pad is made of 600d canvas.
3. Slanted armrests, PU integrated armrest pads.
4. Front wheel rubber tire, aluminum fork.
Rear wheel 20 inches non-inflatable rubber tire, nylon one-piece big wheel.
Large wheels.
5. With rear hand brake function and elbow knuckle brake.
6. Wheelchair back can be folded, armrests can be upturned and fixed, pedals can be upturned.
7. Plastic armrests ring for manual driving of the wheelchair.

Aluminum

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
13/15KG	100KG	91CM	60CM	100CM	45CM	36CM	6 inches	20 inches	/	43CM	91*28*100CM	78*28*70CM

YSM2013



FEATURES:

1. Seat back cushion is made of flame retardant standard oxford nylon fabric mesh cushion, seat back cushion sandwich is 600d canvas.
2. PU one-piece armrest pad.
3. Front wheel 6 inches (15cm)
Rubber tires, high strength steel fork. Rear wheel 20 inches (50cm) solid rubber tire steel wire big wheel.
4. Brake: with rear hand brake function and elbow joint type brake, below the seat surface after braking.
5. Wheelchair back can be folded, armrests can be upturned and fixed, footrests can be upturned. Plastic armrest rings.

Steel

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
12/14KG	100KG	80CM	58CM	90CM	43CM	37CM	6 inches	16 inches	40CM	48CM	78*24*71CM	80*26*73CM

YSM2032



FEATURES:

1. Aluminum alloy frame, lightweight and fast folding.
2. Quick release rear wheel design, up and down in just one second.
3. Latex cushion, comfortable and worm-proof.
4. Adjustable armrests and removable footrests.
5. Three height adjustments of the rear axle.

Quick release design

Height adjust

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
15/17kg	100KG	105CM	70CM	90CM	45CM	42CM	6 inches	12 inches	/	38CM	105*33*73 CM	110*35*76 CM

R

Manual Wheelchair - Standard Wheelchair

YSM2016

TYPE: A



FEATURES:

1. Four brakes with toilet
 2. 1.2mm tube thickness
 3. Optional pedals: iron, aluminum, plastic,
 4. Armrest optional: blow molding, leather,
 5. Cushion optional, leather/checkered/black
- Tires: solid tires

Bedpan

NW: 16KG

TYPE: B



FEATURES:

1. Spraying, aluminum alloy rim
2. One potty, mesh breathable cushion
3. Elevating legs, suitable for orthopedics

Orthopedic leg rests

NW: 18.5KG

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	Folding width	DIM
18.5KG	100KG	103CM	68CM	90CM	46CM	45CM	7 inches	24 inches	45CM	52CM	26CM	94*26*93CM

YSM2017

TYPE: A



FEATURES:

1. Fixed armrest and footrest
2. Comfortable leather seat cushion
3. Foldable and light for easy travel
4. Big rear wheel manual self drive
5. Carbon steel chromed frame

Plating 809

N/G: 19/21KG

TYPE: B



FEATURES:

1. Fixed armrest and footrest
2. Comfortable leather seat cushion
3. Foldable and light for easy travel
4. Big rear wheel manual self drive
5. Carbon steel chromed frame
6. Detachable bedpan
7. Safe rear handle brake

Plating 809

Bedpan

N/G: 18/20KG

TYPE: C



FEATURES:

1. Net backrest sticky pad
2. Black Oxford cloth cushion
3. Gray PU tire plastic bracelet
4. Plastic pedals
5. 8 inch narrow front wheel
6. Plastic guard
7. Soft armrests
8. Single brake

Spraying 809

PU rear tire

N/G: 14/16KG

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
/	100KG	103CM	68CM	85CM	42CM	45CM	8 inches	24 inches	42CM	51CM	103*27*85CM	94*23*89CM

Manual Wheelchair - Standard Wheelchair

YSM2033



FEATURES:

1. Comfortable honeycomb, breathable cushion
2. Portable handbag, foldable & light for easy travel
3. The widened design makes the riding space bigger and more comfortable
4. Carbon steel raw materials reinforcement design

Wider seat

Spraying

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	Folding size	DIM
19/21KG	100KG	94CM	75CM	110CM	56CM	/	8 inches	24 inches	/	/	79*37*94CM	100*26*90CM

YSM2034



FEATURES:

1. The frame is made of steel tube with spraying treatment, round tube frame design, tube diameter 22mm, wall thickness 1.0mm, the frame is beautiful and solid.
2. The handrail can be upturned to be movable and become a triangular structure, which is convenient for users to get on and off the car, and the handrail pad is made of leather sponge material, which is sturdy and durable.
3. Configuration of plastic footrest, unarmed can be turned back and folded to reduce space for easy storage, footrest and adjustable height, suitable for people of different heights, high strength is not easy to break.

Carbon Steel Frame

Spraying

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding width	DIM
19/20KG	100KG	100CM	68CM	85CM	45CM	45CM	7 inches	24 inches	45CM	49CM	26.5CM	93*25*93CM

YSM2046



FEATURES:

1. Spraying the wheelchair
2. 1 pair of detachable front hanging legs
3. Solid PU wheel
4. The armrest can be folded
5. Full mesh breathable cushion
6. Color optional: turquoise, dark blue, black, red

Carbon Steel Frame

MAIN PARAMETERS

N	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
12KG	200KG	85CM	56CM	85CM	45CM	38CM	6 inches	12 inches	41CM	45CM	/	94*28*89CM

YSM2052



FEATURES:

1. Aluminum alloy painted camouflage one-wheel wheelchair
2. Movable folding handle, 20-inch integrated wheel tire
3. Color: camouflage color

Aluminum Frame

Integrated wheel

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
/	/	95CM	60CM	86CM	44CM	37CM	12CM	30CM	/	45CM	81*60*70CM	70*79*26CM

Manual Wheelchair -Reclining Wheelchair

YSM2023



FEATURES:

- 1.Portable commode
- 2.Adjustable armrest and footrest
- 3.Thick and washable bee net
- 4.Hydraulic control backrest angle
- 5.Free rubber solid wheel with manual ring
- 6.Chromed frame, adjustable wheelchair length and height

Hydraulic reclinable

Bedpan

Electroplated

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	92-165CM	68CM	90-120CM	45CM	41CM	7 inches	24 inches	45CM	51CM	92*26*92CM	80*30*70CM

YSM2022



FEATURES:

- 1.Portable commode
- 2.Adjustable armrest and footrest
- 3.Thick and washable bee net
- 4.Hydraulic control backrest angle
- 5.Free rubber solid wheel with manual ring
- 6.Chromed frame, adjustable wheelchair length and height

Hydraulic reclinable

Bedpan

N/G: 20/24KG

rear wheel 24 inches

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	92-165CM	68CM	90-120CM	45CM	41CM	7 inches	/	45CM	51CM	92*26*92CM	80*30*70CM

IYASOCARE®

Manual Wheelchairs - Featured Wheelchairs

YSM2025



FEATURES:

1. PU wheels and hand rims.
2. Adjustable belt handhold ring drive.
3. The user's hands do not come into direct contact with the rear wheel, avoiding the risk of 2 cross-infections.
4. More ergonomic, labor-saving and comfortable, to avoid long-term use caused by frozen shoulder and other health hazards.

Ergonomics

Prevention of cross-infection

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	104CM	68CM	93CM	46CM	40CM	8 inches	13 inches	/	45CM	104*37*76 CM	77*40*78 CM

YSM2059



FEATURES:

1. Spray frame, high back, adjustable and detachable headrest
2. Adjustable front and rear seats, adjustable backrest Angle, adjustable armrest height, removable lifting foot support
3. TPR front wheel, plastic wheel frame inflation quickly remove the rear wheel
4. Rear handle with brake

Car seat

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N	150KG	88CM	64CM	121-133CM	46CM	/	8 inches	24 inches	/	58CM	/	/

YSB01



FEATURES:

1. Foldable frame, front and rear wheels quick detachable installation design
2. Adjustable footrest and armrest design, easy to switch the rest position and get on and off the wheelchair.
3. Full steel frame, sturdy and durable, can carry 150kg.
4. Seat and back cushions are made of glass fibre, breathable and waterproof.
5. The front and rear wheels are balloon wheels, which makes it easy to drive on the beach or wade in the water.

Beach

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	136CM	102CM	94CM	46CM	50CM	12* 7 inches	16* 9 inches	32CM	56CM	/	89*42*84CM

YSB02



FEATURES:

1. Foldable frame design, front and rear wheels quick release design, can be placed in the trunk of the car.
2. Adjustable leg rest, adjustable backrest, maximum tilt 170°.
3. Stainless steel frame, waterproof and durable, with a weight capacity of 100 kg.
4. The cushion is made of breathable and waterproof oxford fabric.
5. Front and rear wheels are balloon wheels, making it easy to drive on the beach or wade in the water.
6. Optional flotation kit for floating in the water.

Beach

Floating

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	123CM	100CM	106CM	50CM	48CM	10 inches	12 inches	48CM	/	99*63*75CM	82.5*45.5*84.5CM

Manual Wheelchairs - Accessories

Wooden tablecloth		UZ Cushion		16-inch Integrated wheel		20-inch low profile tires	
Plastic dining table		Leather headrest		20-inch Integrated wheel		22-inch low profile tires	
Plastic Hard Handrail		Oxford cloth headrest		24-inch Integrated wheel		24-inch low profile tires	
Leather Foam Armrest		A brake		24-inch nylon Integrated wheel			
Plastic Side Guards		Electroplated hand rim ordinary tire wire wheel		20-inch plastic hand-rimmed wire wheel			
Iron side guards		Normal low tire		22-inch plastic hand-rimmed wire wheel			
Plastic foot pedal		Plain Black Gray Tall Tire		24-inch plastic hand-rimmed wire wheel			
Iron footrest		20-inch pu foam tire		20-inch Electroplated Hand Rim Wire Wheel			
Aluminum foot pedal		24-inch pu foam tires		22-inch Electroplated Hand Rim Wire Wheel			
Fork		24-inch pneumatic tires		24-inch Electroplated Hand Rim Wire Wheel			
Square Basin		6-inch rubber black/ pu gray front wheel		16-inch tires (tall tires)			
Chamber pot		7-inch rubber black/ pu gray front wheel		20-inch tall tires			
Leather leg rests		8-inch rubber black/ pu gray front wheel		22-inch tall tires			
Double Mesh Cushion		12-inch rear wheel		24-inch tall tires			

Sport Wheelchair

YSM2101



FEATURES:

- 1.Fixed baffle
- 2.Folding frame
- 3.Quick remove rear wheel
- 4.Oxford fabric seat and back
- 5.Fixed foot pedal,length adjustable
- 6.Light magnesium alloy frame, surface spray

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
14.5KG/19KG	100KG	80CM	64CM	87CM	38CM	38CM	4 inches	24 inches	34CM	43CM	/	70*63*60CM

YSM2103



FEATURES:

- 1.Frame: the combination of welding of high strength aluminum alloy profile forming, fixed structure, safety performance is good, the surface color spray paint processing, beautiful and durable
- 2.Front wheel: 18 inches PU special racing wheel, adopt the whole aluminum alloy brake, safe and reliable
- 3.Rear wheel: 28 inches inflatable wheels, configuration is 36 knife type high strength aluminum circle, inflatable tires, with light weight, high strength aluminum wheel (hand when using direct drive wheelchair)
- 4.Seat cushion is made of 420D oxford cloth stitching forming, soft texture, high tensile strength, the soft cushion can according to the different requirements of different figure athletes on the wheelchair seat height and the thickness of the configuration to accommodate

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
11.8KG/14KG	100 KG	164CM	63CM	68CM	30-34CM	36CM	18 inches	28 inches	/	45CM	/	165*56*70CM

YSM2102



FEATURES:

- 1.Easy folding, after folding can be put in the trunk
- 2.Quick disassemble rear wheel: Press the quick disassemble button and the rear wheel can be disassembled easily
3. Aviation aluminum alloy frame
4. Anti-overturn wheels are provided on both sides of the back of the wheelchair

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
11.6KG/14KG	100KG	85CM	56CM	77CM	36CM	48CM	4 inches	24 inches	/	48CM	85*34*77CM	87*36*79CM

YSM2103C



FEATURES:

- 1.Frame: the combination of welding of high strength aluminum alloy profile forming, fixed structure, safety performance is good, the surface color spray paint processing, beautiful and durable
- 2.Front wheel: 18 inches PU special racing wheel, adopt the whole aluminum alloy brake, safe and reliable
- 3.Rear wheel: 28 inches inflatable wheels, configuration is 36 knife type high strength aluminum circle, inflatable tires, with light weight, high strength aluminum wheel (hand when using direct drive wheelchair)
- 4.Seat cushion is made of 420D oxford cloth stitching forming, soft texture, high tensile strength, the soft cushion can according to the different requirements of different figure athletes on the wheelchair seat height and the thickness of the configuration to accommodate

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
12KG/15KG	100KG	164CM	63CM	68CM	28/30/32 CM	36CM	18 inches	26 inches	/	45CM	/	165*65*70CM

Sport Wheelchair

YSM2104



FEATURES:

1. Easy folding, after folding can be put in the trunk
2. Quick disassemble rear wheel: Press the quick disassemble button and the rear wheel can be disassembled easily
3. Aviation aluminum alloy frame, designed for dance moves
4. Anti-overturn wheels are provided on both sides of the back of the wheelchair

Dance

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	85CM	56CM	75CM	36CM	/	2.4 inches	22 inches	/	50CM	/	85*58*80CM

YSM2106



FEATURES:

1. Easy folding, after folding can be put in the trunk
2. Quick disassemble rear wheel: Press the quick disassemble button and the rear wheel can be disassembled easily
3. Aviation aluminum alloy frame, designed for tennis
4. Anti-overturn wheels are provided on both sides of the back of the wheelchair

Tennis

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	85CM	85CM	68CM	30/32/34/36/38/40CM	/	3.2 inches	26 inches	/	42CM	/	103*51*67CM

YSM2105



FEATURES:

1. Easy folding, after folding can be put in the trunk
2. Quick disassemble rear wheel: Press the quick disassemble button and the rear wheel can be disassembled easily
3. Aviation aluminum alloy frame, designed for basketball guard
4. Anti-overturn wheels are provided on both sides of the back of the wheelchair

Basketball

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	85CM	85CM	64CM	30/32/34/36/38/40CM	/	4 inches	26 inches	/	50CM	/	82*80*86CM

YSM2107



FEATURES:

1. Aluminum wheelchair frame
2. Quick release of rear wheel, convenient storage of wheelchair
3. Wear a seat belt
4. Sturdy plastic footboards

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	75CM	67CM	77CM	36CM	/	4 inches	24 inches	/	48CM	/	76*62*50CM

Sport Wheelchair

YSM2108



FEATURES:

1. Easy folding, after folding can be put in the trunk
2. Quick disassemble rear wheel: Press the quick disassemble button and the rear wheel can be disassembled easily
3. Aviation aluminum alloy frame
4. Anti-overturn wheels are provided on both sides of the back of the wheelchair

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	81.5CM	57-63CM	76.5CM	36/38/40/42CM	40CM	4 inches	24 inches	/	46CM	81.5*35*76.5CM	69*38*85CM

YSM2111



FEATURES:

1. Lightweight design, fashion modeling
2. Carbon fiber frame
3. Anti-overturn wheels are provided on both sides of the back of the wheelchair

Carbon fiber

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	84CM	63-69CM	77CM	36/38/40/42CM	40CM	4 inches	24 inches	/	46CM	/	90*65*78CM

YSM2109



FEATURES:

1. Aluminum foldable frame
2. Fixed armrests, foldable and height-adjustable footrests
3. wider solid front wheels, wider off-road quick release light-weight pneumatic rear wheels
4. Comfortable foam seat and back cushion

All terrain

Offroad

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
12.8KG/15.2KG	100KG	84CM	60CM	77CM	36/38/40/42CM	37CM	4 inches	24 inches	/	48CM	87*33*76 CM	87*36*79 CM

YSM2112



FEATURES:

1. Lightweight design, fashion modeling
2. Aluminum frame
3. Aluminum alloy manual parking brake
4. Cushion fabric selection of high-quality mesh fabric lined with high-quality high resilience sponge

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
16/18KG	100KG	80CM	54CM	76CM	36/38/40/42CM	37CM	5 inches	24 inches	30CM	52CM	/	68*64*58CM

Cerebral Palsy Wheelchair

YSM2201



FEATURES:

1. One-piece Oxford cloth backrest cushion with butterfly seat belt and double brake;
2. ABS guard plate. Feet can be lifted up and removable. Easy to carry;
3. ABS foot pedal, PU handrail, upper and lower three gears adjustable.
4. 135° adjustable back, 90° backrest, adjustable Angle.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	112CM	46CM	100CM	38CM	39CM	6 inches	16 inches	/	46CM	80*25*90CM	83*30*94CM

YSM2206



FEATURES:

1. Aluminum alloy spray frame;
2. 3D mesh surface cushion, PVC wheels;
3. 150° adjustable back, 90° backrest, adjustable Angle.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	75KG	95CM	46CM	112CM	33CM	35CM	6 inches	12 inches	45CM	60-72CM	75*37*90CM	85*51*107CM

YSM2205



FEATURES:

1. Butterfly seat belt and double brake;
2. Feet can be lifted up and removable. Easy to carry;
3. Aluminum alloy front fork, ABS foot pedal, PU handrail, upper and lower three gears adjustable.
4. 150° adjustable back, 90° backrest, adjustable Angle.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	75KG	92CM	42CM	103CM	30CM	30CM	6 inches	12 inches	45CM	49CM	75*37*90CM	80*41*69CM

YSE402

Electric



FEATURES:

1. Aluminum frame lightweight folding design, easy to carry and travel
2. Smart controller, easy to study
3. Detachable footrest
4. Double back handle
5. Flip-up armrest

PARAMETERS:

Climb Angle: 8 degrees
 Motor: 200w*2 brushless motor
 Battery: 24V 10Ah lithium
 Charge time: 4-6 H
 Maximum speed: 6 km/h
 Maximum range: 15 km

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
21/24kg	100KG	100CM	52CM	102CM	39CM	39CM	7 inches	12 inches	28CM	58CM	74*31*69 CM	79*38*75 CM

Children Wheelchair

YSE403

Electric



FEATURES:

1. Aluminum frame lightweight folding design, easy to carry and travel
2. Smart controller, easy to study
3. Detachable footrest
4. Double back handle
5. Flip-up armrest

PARAMETERS:

Climb Angle: 8 degrees
 Motor: 200w*2 worm gear motors
 Battery: 24V 23Ah lithium
 Charge time: 4-6 H
 Maximum speed: 6 km/h
 Maximum range: 20-25 km

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100kg	96CM	56CM	92CM	38CM	40CM	7 inches	12 inches	/	52CM	74*31*69CM	86*58*75CM

YSM2203



FEATURES:

1. The surface is treated with baking paint, which has the characteristics of no fading and strong anti-rust ability.
2. Foldable structure with high safety performance
3. Oxford cloth fabric sewing molding, soft texture, breathable and non-slip
4. backrest armrest can be adjusted, footrest can be detached.

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	100CM	51CM	89CM	35CM	36CM	6 inches	22 inches	41CM	49CM	/	78*41*92CM

YSM2202



FEATURES:

1. Aluminum alloy frame, tube diameter 25MM, wall thickness 1.5MM, baking paint.
2. the footrest height can be adjusted and can be rotated to facilitate easy to get on and off the bike, the whole bike can be folded
3. The seat is high in the front and low in the back, and equipped with a strap-type three-point safety belt to prevent the child from slipping down.
4. 16 inch PU rear wheels, aluminum alloy parking brake.
5. With sunshade, with armrests handrail

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	115CM	57CM	107CM	44CM	35-45CM	8 inches	16 inches	62CM	55-45CM	120*38*35CM	120*38*35CM

YSM2204



FEATURES:

1. An adaptable aluminium manual wheelchair designed for children with disabilities.
2. The seat width and seat depth can be adjusted. With the growth of children, the seat width can be adjusted to 25cm, 30cm and 35cm
3. Seat depth can be adjusted to 25cm, 30cm, 35cm three kinds of seat depth.

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	75KG	92CM	60CM	92CM	25-35CM	25-35CM	4 inches	20 inches	/	36-41CM	92*36*87CM	78*41*92CM

Children Wheelchair

YSM2208



FEATURES:

1. Comfortable and breathable cushion
2. Aluminium alloy frame
3. Solid rubber wheels with manual ring
4. Rear brake handle
5. Removable footrest

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	75KG	75CM	48CM	94CM	30CM	30CM	4 inches	20 inches	38CM	46CM	65*30*69CM	71*33*69CM

YSM-PONY



FEATURES:

1. the frame aluminum alloy material, the surface is painted pearl white, the whole vehicle can be folded.
2. backrest 3 gear can be adjusted back 90-110 degrees, 5 inch universal front wheel, 5 inch directional rear wheel, treadle brake.
3. 3 gear sitting depth can be adjusted.
4. foot support height is adjustable.
5. the head pillow height is adjustable.
6. hips with 3-point seat belt.
7. maximum load 75KG.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	75KG	96CM	40CM	94CM	33CM	25CM	/ inches	/ inches	/	48CM	103*40*34CM	110*47*42CM

YSM2209



FEATURES:

1. Comfortable and breathable cushion
2. Aluminium alloy frame
3. Solid rubber wheels with manual ring
4. Adjustable handrest
5. Removable footrest

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	75KG	95CM	46CM	90CM	35CM	36CM	6/ inches	22/ inches	38CM	48CM	83*28*90CM	85*29*92CM

IYASOCARE®

Showers And Commode Wheelchair

YSM2007



FEATURES:

- 1, Three-in-one wheelchair, commode chair, shower chair, manual wheelchair
- 2, Detachable commode for easy cleaning
- 3, Detachable footrest, easy to fold
- 4, Adjustable leather armrests
- 5, Four small wheels with independent brakes

YSM2057



FEATURES:

1. Aluminum alloy frame with folding handle
2. Fabric material: luxurious sofa cushion
3. Wheel Material: Shock Absorbing PU
4. Waterproofing

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
13/15KG	100KG	100CM	54CM	93CM	41CM	41CM	/	/	/	47CM	81*26*86CM	67*60*66CM

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
13/15KG	100KG	100CM	54CM	93CM	41CM	41CM	/	/	/	47CM	81*26*86CM	67*60*66CM

YSM2058



FEATURES:

1. Aluminum alloy frame with folding handle
2. Fabric material: luxurious sofa cushion
3. Wheel Material: Shock Absorbing PU
4. Waterproofing

MAIN PARAMETERS

NW	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
13KG	100 KG	91CM	57CM	90CM	/	/	6 inches	20 inches	40CM	/	91*25*70CM	71*29*70CM

IYASOCARE®

Transfer Wheelchair

YSM2301



FEATURES:

1. Medical mute wheel, hub brake
2. Waterproof design, easy to shower and cool
3. Equipped with infusion stand jack, can be disassembled and replaced freely

Manual

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	59CM	54CM	91-106CM	49CM	45CM	12 inches	7.5 inches	37CM	45-62CM	/	36*61*74CM

YSM2304



FEATURES:

1. Transfer lift chair
2. Applicable area: home and hospitable use
3. Transfer patients and disabled or elder people, waterproof can be use for taking shower
4. Hydraulic power height lift adjustable

Hydraulic

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100 KG	74CM	59CM	91-106CM	49CM	40CM	7 cm	5 cm	46CM	45-66CM	/	99*70*40CM

YSM2302



FEATURES:

1. Transfer lift chair.
2. Applicable area: home and hospitable use.
3. Transfer patients and disabled or elder people, waterproof can be use for taking shower.
4. Electric power height adjustable.

Electric

PARAMETERS:

Color: white and black
 Material: metal, Plastic, Aluminum Alloy/Foamed plastics
 Drive mode: Motor driven by Lithium battery
 Input voltage 110-220v
 Output voltage 24v

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	/	/	/	51CM	/	/	/	/	/	/	86*65*30CM

YSM2305



FEATURES:

1. Transfer lift chair
2. Applicable area: home and hospitable use
3. Transfer patients and disabled or elder people, waterproof can be use for taking shower
4. Hydraulic power height lift adjustable

Hydraulic

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	100KG	85CM	54CM	92-112CM	47CM	39CM	/	/	46CM	45-65CM	/	73*62*39CM

Walking Frame

YSW7101



FEATURES:

1. Standing exercise walking aid
2. Function: standing exercise/ suitable for home use/flexible moving
3. Accessories that can be added
The size can be customized according to customer needs
4. Paraplegic patients can easily stand up, exercise easily, and move freely after standing with the help of double-arm tension and a controllable hydraulic power assist system. Improve blood circulation, prevent complications such as muscle atrophy, bedsores, and improve the quality of life of patients.

PARAMETERS:

Armrests height: 103cm(high armrest)/73cm(low armrest)
 Armrests length: 30cm(high armrest)/27cm(low armrest)
 Inner width of upper armrest: 46cm

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	100CM	65CM	140CM	23CM	45CM	10 inches	3 inches	/	50CM	/	112*68*65CM

YSW7102C



FEATURES:

1. 6-wheels stability without overturning
2. 10-level height adjustment is suitable for different groups of people
3. Improve blood circulation, prevent complications such as muscle atrophy, bedsores, and improve the quality of life of patients.

PARAMETERS:

Height adjustment: 70-90cm (10 gears)

Children size

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
9/11KG	75KG	90CM	40CM	75CM	/	/	/	/	/	/	/	48*69*25CM

YSW7102



FEATURES:

1. 6-wheels stability without overturning
2. 10-level height adjustment is suitable for different groups of people
3. Improve blood circulation, prevent complications such as muscle atrophy, bedsores, and improve the quality of life of patients.

PARAMETERS:

Height adjustment: 116-135cm (10 gears)

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
14.3KG/16.5KG	100KG	108CM	72CM	135CM	/	/	4.2 inches	4.2 inches	/	/	/	75*25*89CM

YSW7103



FEATURES:

1. Standing exercise walking aid
2. Function: standing exercise/ suitable for home use/flexible moving
3. Accessories that can be added
The size can be customized according to customer needs
4. Paraplegic patients can easily stand up, exercise easily, and move freely after standing with the help of double-arm tension and a controllable hydraulic power assist system. Improve blood circulation, prevent complications such as muscle atrophy, bedsores, and improve the quality of life of patients.

PARAMETERS:

Armrests height: 103cm(high armrest)/
 73cm(low armrest)
 Armrests length: 30cm(high armrest)/
 27cm(low armrest)
 Inner width of upper armrest: 46cm

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
65KG	150KG	90CM	65CM	130CM	23CM	45CM	10 inches	3 inches	/	50CM	/	112*68*65CM

Rollator

YSW7005



FEATURES:

1. PU material for front and rear wheel
2. Self-locking handle brake
3. Double folding frame

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
7.4kg	100KG	68CM	58CM	84-93CM	46CM	32*35CM	/	/	/	/	/	88*60*25CM

YSW7010



FEATURES:

1. PU material for front and rear wheel
2. Self-locking handle brake
3. Double folding frame

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
9.6kg	130KG	95CM	61CM	85-97CM	45CM	/	/	/	/	54CM	70*28*85CM	73*30*88CM

YSW7009



FEATURES:

1. PU material for front and rear wheel
2. Self-locking handle brake
3. Double folding frame

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
7.4KG	136KG	67CM	55.5CM	83.5-96CM	45CM	/	/	/	/	54CM	/	42*23*69CM

YSW7011

Overspeed protection



FEATURES:

1. PU material for front and rear wheel
2. Self-locking handle brake
3. Double folding frame
4. Resistance adjustment for safe speed-limited travel without electric power

Resistance is generated by a magnetic field to create resistance without the need for electric power. This technology can limit the speed of movement of the walker, provide stable support for the user, and assist the user to walk within the controllable speed range, and the faster the walking speed, the greater the resistance.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
/	130KG	74CM	65CM	79-87CM	45CM	30CM	8	8	/	52CM	74*32*85CM	75*34*90CM

Rollator

YSW7013



FEATURES:

- 1, The front and rear wheels are solid universal wheels, easy to cope with different road surfaces
- 2, Three height adjustments of the armrests to meet the use of people of different heights.
- 3, Quick foldable design, easy to carry and store.
- 4, Triangular frame design, stable and anti-tip over.
- 5, Breathable cushion and sponge backrest make the ride more comfortable.

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
8/9kg	120KG	47CM	51CM	89-96CM	43CM	/	7 inches	7 inches	/	/	74*32*85CM	76*53*35CM

YSW7015

Overspeed protection



FEATURES:

1. No batteries, four resistance speed adjustment of the front wheel, suitable for different use of the field. Scene
2. Aluminum alloy frame, easy to fold and portable.
3. Armrest height adjustable, suitable for different people.
4. Independent brakes on front and rear wheels, safer parking.
5. Removable seat plate, ready to use, ready to dismantle.

MAIN PARAMETERS

NW	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
17.6kg	100KG	69CM	65CM	83-96CM	44CM	11CM	2.5 inches	5.5 inches	/	54CM	25*21*83CM	28*23*85CM

YSW7014

Electric



FEATURES:

1. Power wheelchair, walker, rollator 3 in 1
2. Aluminum alloy frame can be folded easily and quickly.
3. Ternary lithium battery and brushless motor, lighter weight and longer range.
4. Nylon cushion, comfortable and breathable.
5. Height adjustable backrest, suitable for different heights.
6. Five speed control, so that the elderly can also be easily operated.

PARAMETERS:

Charge time: 2-3H
 Driving Range: 10-20KM
 Max Speed: 6KM/h
 Climbing ability: 15°
 Motor: 250W*2 brushless motor
 Battery: 24V 5 /10Ah ternary lithium-ion battery
 Parking Brake: Manual brake and electronic parking brake
 Cushion: Black

MAIN PARAMETERS

NW	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
17.6KG	100KG	100CM	66CM	80-90CM	45CM	25CM	8 inches	8 inches	/	/	77*29*90CM	73*34*81CM

YSW7016

Electric



FEATURES:

1. Power wheelchair, walker, rollator 3 in 1
2. Aluminum alloy frame can be folded easily and quickly.
3. Ternary lithium battery and brushless motor, lighter weight and longer range.
4. Nylon cushion, comfortable and breathable.
5. Height adjustable backrest, suitable for different heights.
6. Five forward gears, speed: 2KM/H, 3KM/H, 4KM/H, 5KM/H, 6KM/H. One reverse gear: 1.0~2.2KM/H

PARAMETERS:

Charge time: 5-6H
 Driving Range: 10-20KM
 Max Speed: 6KM/h
 Climbing ability: 6°
 Motor: 250W*2 brushless motor
 Battery: 24V 5 /10Ah ternary lithium-ion battery
 Parking Brake: Manual brake and electronic parking brake
 Cushion: Black

MAIN PARAMETERS

N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
19/24.7KG	115KG	100CM	70CM	76-96CM	45CM	35CM	10 inches	8 inches	/	/	73.5*32.5*76CM	74*36*85CM

Walking Aid

YSW7002



FEATURES:

1. Made of steel, it will not scratch the surface and never rust;
2. The height of the legs is adjustable in 8 gears and 8CM, which is suitable for different heights (there is a slight error due to different production batches);
3. The handle is made of high-density PVC material, which is non-slip and does not break;
4. Enlarged and thickened natural rubber foot pads are more durable and non-slip, and built-in metal gaskets enhance their durability;
5. When not in use, it is very convenient to fold and occupies a small area.

Instructions

1. Fold and open, and the middle marble pops out to fix;
2. Press the 4-leg marbles in turn, and adjust the 4 legs to the same appropriate height;
3. Hold the handle tightly and move forward alternately with the left hand and the right hand.

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	52CM	61CM	77-95CM	/	/	/	/	/	/	80*61*10CM	54*38*78CM

YSW7004



FEATURES:

1. Fold and open, and the middle marble pops out to fix;
2. Press the 4-leg marbles in turn, and adjust the 4 legs to the same appropriate height;
3. Hold the handle tightly and move forward alternately with the left hand and the right hand.

PARAMETERS:

- 1.8-level height adjustment,
- 2.Two 13cm diameter front casters,
- 3.Aluminum alloy oxidation treatment,
- 4.Two rear anti-slip feet with built-in metal spacers.
- 5.H-Link Crossbar, Marble Folding

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	52CM	61CM	77-95CM	/	/	/	/	/	/	80*61*10CM	54*38*78CM

YSW7003



FEATURES:

- 1.8-level height adjustment,
- 2.Two 13cm diameter front casters,
- 3.Aluminum alloy oxidation treatment,
- 4.Two rear anti-slip feet with built-in metal spacers.
- 5.H-Link Crossbar, Marble Folding

Instructions:

1. Fold and open, and the middle marble pops out to fix;
2. Press the 4-leg marbles in turn, and adjust the 4 legs to the same appropriate height;
3. Hold the handle tightly and move forward alternately with the left hand and the right hand.

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	51CM	61CM	77-94CM	/	/	/	/	/	/	80*61*10CM	54*38*78 CM

YSW7008



FEATURES:

- 1.5-level height adjustment,
- 2.Stable front casters,
- 3.Steel framed structure,
- 4.Two rear anti-slip feet ,
- 5.X frame for easy folding.

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	150KG	56CM	61CM	81-90CM	43CM	/	/	/	/	49-58CM	84*61*26CM	54*38*78CM

Crutch And Cane



YSC6001

height: 73-96cm
N.W:0.9KG



YSC6002

height: 64-87cm
N.W:0.9KG



YSC6003

height: 70-93cm
N.W:0.35KG



YSC6005

height: 102.5-144.5cm
N.W:1KG
Loading: 200KG



YSC6007

height: 88-109.5cm
NW: 0.88KG



YSC6008

Product size: 27* (84-92) cm
Adjustable height: 84-92cm
Seat deep: 24cm
Seat height: 47-56cm
NW: 1.08KG
Loading: 150KG
Material: Aluminium



YSC6009

Product height: 84-142cm
Applicable user height:
140-180cm
NW: 1.05KG
Loading: 130KG
Material: Aluminium



YSC6010

Product height: 95-118cm
Folding size: 54*19*12.5cm
Applicable user height:
150-180cm
NW: 0.54KG
Loading: 136KG
Material: Aluminium



YSC6011

Applicable user height: 150-180cm
NW: 0.6KG
Loading: 130KG
Material: Aluminium
Carton size: 82*36*36cm
10 pairs/ctn
N.W/G.W: 12/13.1kg



YSC6018

Cruth model height: 84cm
Seat model size: 26*35*72cm
Seat height: 50cm
NW: 0.9KG
Loading: 100KG
Material: Aluminium frame
and ABS panel

Shower And Commode Chair

Y>>> YSB8001



MAIN PARAMETERS												
N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
15/16 kg	/	/	/	69~81.5CM	50 CM	42CM	/	/	/	39~51.5CM	/	65*52*42CM

Y>>> YSB8002



MAIN PARAMETERS												
N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
3/7KG	/	/	/	71~83.5CM	47CM	46CM	/	/	/	41~53.5CM	/	51*38*29CM

Y>>> YSB8004



MAIN PARAMETERS												
N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
6/13KG	/	/	/	78~88CM	44CM	52CM	/	/	/	45~55CM	/	52*34*75CM

Y>>> YSB8005



MAIN PARAMETERS												
N/G	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
3.4/4.4KG	/	70CM	57CM	84CM	50CM	50CM	/	/	/	/	/	52*34*75CM

Shower And Commode Chair

YYSB8008



YYSB8010



MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	4.2KG	/	46CM	48CM	66-76CM	46CM	/	/	/	35-45CM	/	87*58*61CM

MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	4.5KG	/	/	/	64-74CM	43CM	43CM	/	/	36-46M	/	52*34*75CM

YYSB8011



MAIN PARAMETERS

NG	Load	overall length	overall width	overall height	seat width	seat depth	front wheel	rear wheel	backrest height	seat height	folding size	DIM
N/G	3.3KG	/	48CM	78CM	89CM	43CM	/	/	/	45CM	78*53*17CM	82*58*56CM

IYASOCARE®

Support Backrest

YSK1



FEATURES:

- 1, Removable pillow: easily removable, easy to clean, durable and comfortable.
- 2, 6 adjustable gears, angle 0-80° can be adjusted according to the comfort level of the backrest tilt .
- 3, Foldable non-slip armrests, the backrest frame is equipped with armrest support, which can make the user feel comfortable, and the armrests are parallel to the backrest, which will not pinch the user.
- 4, Foldable when not in use, taking up less space and more convenient.
- 5, Breathable waterproof cloth, easy to clean

MAIN PARAMETERS

N/G	Load	overall length	overall width	Frame material	Maximum angle	Motor power	Cushion Size	DIM
2.1/2.55KG	100KG	68CM	58CM	STEEL	80°	/	37*37CM	66*7*70CM

YSK2



FEATURES:

- 1, Removable pillow: easily removable, easy to clean, durable and comfortable.
- 2, 6 adjustable gears, angle 0-80° can be adjusted according to the comfort level of the backrest tilt .
- 3, Foldable non-slip armrests, the backrest frame is equipped with armrest support, which can make the user feel comfortable, and the armrests are parallel to the backrest, which will not pinch the user.
- 4, Foldable when not in use, taking up less space and more convenient.
- 5, Breathable waterproof cloth, easy to clean

MAIN PARAMETERS

N/G	Load	overall length	overall width	Frame material	Maximum angle	Motor power	Cushion Size	DIM
2.1/2.55KG	100KG	68CM	58CM	STEEL	80°	/	37*37CM	66*7*70CM

YSK3



FEATURES:

- 1, Removable pillow: easily removable, easy to clean, durable and comfortable.
- 2, 6 adjustable gears, angle 2-80° can be adjusted according to the comfort level of the backrest tilt .
- 3, Foldable non-slip armrests, the backrest frame is equipped with armrest support, which can make the user feel comfortable, and the armrests are parallel to the backrest, which will not pinch the user.
- 4, Foldable when not in use, taking up less space and more convenient.
- 5, Breathable waterproof cloth, easy to clean

MAIN PARAMETERS

N/G	Load	overall length	overall width	Frame material	Maximum angle	Motor power	Cushion Size	DIM
2.1/2.55KG	100KG	68CM	58CM	STEEL	80°	/	37*37CM	66*7*70CM

YSK4



FEATURES:

- 1, Detachable Pillow: Detachable, easy to clean, durable and comfortable.
- 2, Electric adjustable, angle 0-65° can be adjusted according to the comfort of the backrest tilt.
- 3, Foldable non-slip armrests, backrest frame with armrest support, can make users feel comfortable, and armrests and backrest parallel, will not pinch the user.
Linkage
- 4, It can be folded when not in use, occupying less space and more convenient.

Electric

MAIN PARAMETERS

N/G	Load	overall length	overall width	Frame material	Maximum angle	Motor power	Cushion Size	DIM
14/16KG	100KG	190CM	75CM	STEEL	65°	50W	/	80*82*16CM